

## SOFT PLUS ANESTHESIA BREATHING MASKS



Item Code: MDX 005445

Brand: ICU Medical

### **Description:**

Anesthesia Breathing Mask, Medium Adult, Inflation Valve, Soft Cush, Dark Blue Hook Ring, 30/cs

Specification Sheet  
Item No: 10

## VYAIRE MEDICAL VITAL SIGNS ANESTHESIA CIRCUITS



Item Code: CRF G96LXXXA

Brand: Vyair Medical, Inc.

### **Description:**

Anesthesia Circuit, Pediatric, 108" Expandable Hose, 1L Breathing Bag,  
15/cs

Specification Sheet  
Item No: 17

## VYAIRE MEDICAL ANESTHESIA BREATHING BAGS



Item Code: GEM M1036387


Brand: Vyaire Medical, Inc.

### **Description:**

2L Bag, Disposable, 20/pk

Veterinary anesthesia breathing bags are a critical component to your anesthesia machine. Breathing bags provide a visual monitor of the patient's breathing, protects the patient from excessive pressure and hold a collection of respiratory gases used during patient's spontaneous breathing.

## Dry Blood Collection Bag 300 ml



Item Code: J0520B

Brand: JorVet

### Description:

A wide range of Blood Transfusion Products allows for all species to be recipients of blood transfusion medicine. Canines tolerate standard human transfusion products. Feline patients require much smaller 50-150ml bags. Equine patients need 1-4 liter bags to do the job adequately. Blood component therapy is an expanding field. Plasma and concentrated RBCs can be separated, and specific therapeutic modalities can be implemented. The addition of transfer bags and dry collection bags make this an easy procedure.

Dry Blood Collection Bag. Provide your own anticoagulant for the use of these blood bags. Two spike ports with donor tubing with "Y" injection port. Includes directions for use.

Collection Bag. 300ml.

*Actual Weight: .35 lbs.*



## JorVet Soda Lime, 3 Lb.



Item Code: J0553

Brand: Jorvet

### Specification:

JorVet Soda Lime. Carbon dioxide absorbent granules. JorVet soda lime is a high-grade granular product that effectively absorbs carbon dioxide that is produced as an anesthesia by-product. It is a mixture of calcium hydroxide  $\text{Ca(OH)}_2$  also known as lime plus a smaller amount of sodium hydroxide  $\text{NaOH}$ . Upon exposure to carbon dioxide  $\text{CO}_2$ , the calcium hydroxide converts to calcium carbonate  $\text{CaCO}_3$ . This prevents  $\text{CO}_2$  buildup in a re-breathing circuit. The sodium hydroxide acts as a catalyst in the process. No potassium hydroxide ( $\text{KOH}$ ) is present in JorVet soda lime; this helps to avoid the formation of toxic compounds such as Compound A and carbon monoxide. JorVet soda lime has a pH color change indicator that changes from white to violet when the soda lime needs to be changed. It will eventually return to its original color if not changed. Therefore, a good rule of thumb is to replace the soda lime every time the vaporizer is filled.

3lb (1.36kg) Foil lined bag.

# Bag, Waste, Dog 200PK



Item Code: TCO28600

Brand: TATCO

## Specification:

Dog waste bags on a roll help keep your areas clean and free from dog waste. Opaque design hides containment inside for added privacy. Durable construction prevents tearing for sanitary use. Universal design allows convenient use with virtually all waste station bag dispensers for rolled bags. Waste bags are perfect for use in parks, residential complexes, office landscapes and more.

## Biopsy Punch, 2.5, 25/bx



Item Code: HCP BP25

Brand: EDM3 Company LLC

### Specification:

- One-piece seamless stainless-steel blade.
- Hollow plastic sure grip handle for easy sample extraction.
- Safety capped for your protection.
- Color coded by size.
- Sterility monitors on every punch.
- See-through, tamper evident sterile packaging.

# 18OZ Blade Wash



Item Code: 76300103005

Brand: NEWELL BRANDS DISTRIBUTION

## **Specification:**

18 OZ, Blade Wash, Flushes Away Hair Buildup Trapped Between Upper Lower Cutting Blades, Ideal for Removing Factory Applied Preservative from New Blades, Always Apply Oster Blade Lube Premium Clipper Oil After Using Blade Wash, Spray

## Pocket Pet Carriers for Small Animals, Disposable, Cardboard, 12 Per Case



Item Code: PET-PC

Brand: OASIS

### Specification:

UPC: 647367885859200 lb test white cardboard 10 in. x 6.75 in. x 5.5 in. Safe, no staples, strong glued joints and extra ventilation. Recyclable, approximately 70% post-consumer recycled product. Easy assembly instructions on carriers. Easy one-piece construction. Sold by the dozen.

# Welborn Cat Catcher



Item Code: J0203

Brand: Jorvet

## **Specification:**

The Welbourn cat catcher is designed to aid the veterinarian in safe and easy handling of even the most unruly animals. The Welbourn cat catcher is constructed of tough nylon webbing and can be opened or closed quickly as needed.

Welbourn Cat Catcher.

Specification Sheet  
Item No: 40

Cast Padding Undercast Sof-Rol® 2 Inch X 4 Yard Rayon NonSterile



Item Description: Sof-Rol® Cast Padding Undercast, 2" x 4yd, Cotton, Sterile, 20 rl/cs  
Supplier Item Code: 9052S  
NDC Item Code: BSN 9052S  
Unit of Measure: CS

Product Specifications

Manufacturer #	9052S
Brand	Sof-Rol®
Manufacturer	BSN Medical
Country of Origin	United Kingdom
Application	Cast Padding
Casting Compatibility	Natural Casting Materials / Synthetic
Color	White
Dimensions	2 Inch X 4 Yard
HCPCS	A6441
Material	Rayon
Sterility	NonSterile
Type	Undercast

Specification Sheet  
Item No: 41

Swabstick Dukal™ Cotton Tip Wood Shaft 6 Inch NonSterile 100 per Pack



The applicators features a soft cotton tip that is absorbent and won't lint. Made of a high quality hard wood, these applicators offer superior performance. The sterile 2's are available in an easy to peel package, and the non-sterile are packaged in bags of 100 with a sterilization indicator strip on the outside.

Item Description: Applicator, 6" Cotton Tip, Non-Sterile, 100/bg, 10  
bg/bx, 10 bx/cs (36 cs/plt)

Supplier Item Code: 9006

NDC Item Code: DUK 9006

Unit of Measure: CS

Product Specifications

Manufacturer #	9006
Brand	Dukal™
Manufacturer	Dukal
Country of Origin	Unknown
Application	Swabstick
FSA Eligible - Primary UOM	Yes
Length	6 Inch
Number per Pack	100 per Pack
Shaft Material	Wood Shaft
Sterility	NonSterile
Tip Material	Cotton Tip



Zep Dynamic A II Aerosol Disinfectant, 16 oz. Aerosol Spray, 12 Cans/Case - 351501



Product Description

Zepdynamic A II inactivates HIV-1 (AIDS) virus and Herpes Simplex Types 1 and 2 viruses on precleaned environmental surfaces or objects previously soiled with blood or bodily fluids. Kills TB, Salmonella, Paratyphoid B, Streptococcus, canine-related and other viruses. Kills Athlete`s Foot fungus, Escherichia coli, and Staphylococcus aureus. Controls mold and mildew. Dries instantly leaving no oily film. EPA registered.

Specifications

Size	16 oz
------	-------

Type	Spray Disinfectant
Disinfectant	Quaternary Ammonia
Scent	Citrus
Size Type	16 oz. Aerosol Can
Applicable Materials	Ceramics, Enamel, Fiberglass, Grout, Laminates, Metal, Plastic, Vinyl

Helps Remove	Bacteria, Dust, Fungus, Germs, Grease, Grime, Mildew, Mold, Viruses
Description	Dynamic A II Aerosol Disinfectant
Package Quantity	12 Cans
Manufacturers Part Number	351501
Brand	Zep

## Penrose Surgical Wound Drain, 12" x 5/8", Latex, Sterile, 10/bx, 10 bx/cs



**Brand: Cardinal Health**

**Item no: CHC 30414-058**

### Specifications:

<b>Manufacturer</b>	Cardinal
<b>Country of Origin</b>	TAA Compliant
<b>Application</b>	Wound Drain Tube
<b>Attachments</b>	Safety Pin
<b>Length</b>	18 Inch Length
<b>Material</b>	Latex
<b>Size</b>	5/8 Inch
<b>Sterility</b>	Flat Type
<b>Style</b>	Penrose Style
<b>Type</b>	Flat Type
<b>UNSPSC Code</b>	42295453
<b>X-Ray Compatibility</b>	Consist of a flexible rubber tube which is left inside the site of the surgery to promote drainage of blood, lymph, and other fluids to help reduce the risk of infection

## Drape Sheet, Non-Sterile, Non-Fenestrated, Bulk Packed, 300/cs



**Brand: Busse Hospital Disposables, Inc.**

**Item no: BUS 3270**

Manufacturer #	3270
Manufacturer	Busse Hospital Disposables
Country of Origin	United States
Application	General Purpose Drape
Color	Blue
Drape Name	Embossed Drape
Extremity Cover	Without Covers
Fluid Collection Type	Without Pouch
Instrument Holders	Without Instrument Holder
Reinforcement Type	Without Reinforcement
Securing Method	Without Adhesive
Sterility	NonSterile
Surgical Area	Nonfenestrated
Tube and Cord Holders	Without Cord and Tube Holders
UNSPSC Code	42131701

## Epsom Salt, 16 oz




**Brand: Cumberland Swan/Vi-Jon, Inc.**

**Item no: CUM 1000036823**

### **Description:**

Cumberland Swan/Vi-Jon 1000036823, CUMBERLAND SWAN EPSOM SALTS Epsom Salt, 16 oz, UPC#3105102, 12/cs

# Fecal Loop



## Description

Inexpensive plastic loop with slotted ends for quick and easy fecal samples.

Reusable.

A small, pencil-sized, flexible rod with a loop on each end. This can easily be inserted into the rectum and a small amount of fecal material trapped in the loop.

- Clean, disinfect, and use repeatedly
- Two convenient sizes available

# Cat and Wildlife Glove, Medium



Top Grain Leather Lined with Spectra

Kelvar Material Double Layered

## Description

Designed for maximum dexterity and feel (often missing in protective gloves). Provides a very high level of penetration protection through a combination of special materials. The hand features high quality top grain elk hide completely lined with new Spectra, which is ten times stronger than Kevlar. The gauntlet's double layer of Dupont Kevlar Knit resists sharp teeth, talons and claws, covering up past the elbow. In addition to penetration resistance, the outstanding feature is you hardly know you have them on. Perfect for handling creatures that can not inflict crushing type bites. Excellent for cats, birds, and small animals. These gloves are NOT bite proof. Please use with caution.

Medium

# ENCORE<sup>®</sup> Sensi-Touch<sup>®</sup> PF

## S U R G I C A L   G L O V E S

### Latex powder-free surgical gloves with a strong formulation, for everyday hand protection

- **Secure personal protection:** ENCORE<sup>®</sup> Sensi-Touch<sup>®</sup> PF surgical gloves are powder-free and feature a strong glove formulation, granting everyday protection you can depend on for general surgical procedures
- **Reliable comfort and fit:** These sterile surgical gloves have been designed using Ansell's proprietary anatomic formers to offer comfort and dexterous hand movement
- **Effortless donning:** Thanks to their polymer-coated internal finish, these dependable surgical gloves are both fast and easy to don
- **Reliable grip:** These latex surgical gloves also feature a smooth external surface, with a micro-textured finish, for secure instrument handling
- **Chemo-tested:** They have been tested for use with chemo-drugs, in accordance with ASTM D6978 requirements

**CAUTION:** Safe use of these gloves by latex-sensitized individuals has not been established – this product contains natural rubber latex, which may cause allergic reactions

### Recommended For

- Brief, low-risk examination procedures
- Sterile procedures in the O.R.
- Use as an outer glove when double gloving
- Sterile procedures outside of the O.R.



### ENCORE<sup>®</sup> The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.

# ENCORE<sup>®</sup> Sensi-Touch<sup>®</sup> PF

## SURGICAL GLOVES

### PRODUCT INFORMATION

Material	Natural Rubber Latex
Powder Content	Powder-Free
Color	Natural
Thickness	Standard
Cuff Style	Beaded
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated
Allergy Prevention	None
Grip Level	Moderate
Double Gloving Recommendation	Outer glove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	ASTM D3577, REACH, EN 420:2003 + A1:2009
Quality / Environmental Standards	ISO 13485, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Diethyldithiocarbamate (ZDEC), Zinc Dibutyldithiocarbamate (ZDBC)

### ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	7821PF	7822PF	7823PF	7824PF	7825PF	7826PF	7827PF	7828PF

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
 Sandel US Tel: 866 764 3327  
 Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: [csiac@ansell.com](mailto:csiac@ansell.com)  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
 Tel: +91 22 38124500  
 Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
 Tel: +86(21) 3827 5005  
 Email: [infochina@ansell.com](mailto:infochina@ansell.com)



# GAMMEX® Non-Latex

## SURGICAL GLOVES

### Neoprene, chemical accelerator-free surgical gloves offer advanced allergy protection

- **Advanced allergy protection:** GAMMEX® Non-Latex neoprene surgical gloves are synthetic and eliminate risks of Type I latex allergy. These sterile gloves are made 100% free of chemical accelerators and minimize any risk of sensitization to allergens leading to allergic contact dermatitis.
- **Proven durability:** Strong neoprene glove film deliver the added hand protection needed for the most rigorous procedures
- **Confident grip:** These sterile gloves feature an external micro-textured finish for firmer, more confident instrument handling
- **Effortless donning:** Ansell's donning aid technology includes an internal polymer coating that facilitate easy dry or damp donning
- **Chemo-tested:** In addition, these long-cuff gloves are tested against ASTM D6978 and recommended for the safe handling of chemotherapy drugs\*



\*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

### Recommended For

- Damp donning
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients
- Tested for use with chemotherapy drugs

### FEATURED TECHNOLOGIES



### GAMMEX® Innovation that Protects

Driven by innovation, GAMMEX® surgical gloves are engineered to respond to the evolving user needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomic comfort and superior protection, healthcare workers around the world entrust their hands to ours.

# GAMMEX<sup>®</sup> Non-Latex

## SURGICAL GLOVES

### PRODUCT INFORMATION

Material	Neoprene (Polychloroprene)
Powder Content	Powder-Free
Color	Green
Thickness	Standard
Cuff Style	Straight with SUREFIT <sup>™</sup> Technology
External Glove Surface	Smooth with Micro-Textured Finish
Internal Glove Surface	Polymer Coated with DERMASHIELD <sup>™</sup> Technology
Allergy Prevention	Latex (Type I), Chemical (Type IV)
Grip Level	Moderate
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
Freedom from Holes (Inspection level I)	0.65 AQL
Sterile	Yes
Product Standards	AS/NZS 4179, ASTM D3577, ASTM D6978, ASTM D7160, ASTM F1670, ASTM F739, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282, JIS T9107
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001, Korean GMP
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	None

### ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	340006055	340006060	340006065	340006070	340006075	340006080	340006085	340006090

Contact your Ansell representative for ordering or more information.

**North America**  
 US Tel: 800 952 9916  
 CA Tel: 844 494 7854  
 Email: [insidesalesus@ansell.com](mailto:insidesalesus@ansell.com)  
 Sandel US Tel: 866 764 3327  
 Sandel email: [sandel-cs@ansell.com](mailto:sandel-cs@ansell.com)

**Central & South America**  
 Tel: +52(442) 296 20 50  
 Email: [csiac@ansell.com](mailto:csiac@ansell.com)  
**Brazil**  
 Tel: +55-11-3356-3100  
 Email: [latam.medical@ansell.com](mailto:latam.medical@ansell.com)

**Europe, Middle East & Africa**  
 Tel: +32 (0) 2 528 74 00  
 Email: [info@ansell.eu](mailto:info@ansell.eu)  
**Asia Pacific**  
 Tel: +603 8310 6688  
 Email: [apac.medical@ansell.com](mailto:apac.medical@ansell.com)

**India**  
 Tel: +91 22 38124500  
 Email: [info.india@ansell.com](mailto:info.india@ansell.com)  
**Australia & New Zealand**  
 Tel: +61 3 9270 7270  
 Email: [protection@ap.ansell.com](mailto:protection@ap.ansell.com)  
**China**  
 Tel: +86(21) 3827 5005  
 Email: [infochina@ansell.com](mailto:infochina@ansell.com)

## Powder-Free Latex Exam Gloves

Medgluv Latex Powder-Free Exam Gloves have excellent barrier properties combined with great fit, feel and are powder-free to minimize skin irritation. They offer a soft and comfortable feel with outstanding tactile sensitivity and dexterity.

### Description:

- This 6.5 mil textured glove ensures a secure grip while the superior elasticity provides unequal comfort, fit, strength and tactile sensitivity
- Standard Malaysian Gloves (SMG) certification assurance for quality low-protein latex exam glove.
- Offer excellent sensitivity touch perception and dependable protection
- Minimal low-protein level latex helps reduce latex sensitivity to wearers.
- Biodegradable and Environment friendly.

Item	Description	Size	Packaging
MG100XS	Medgluv Latex Exam Glove	X-Small	100 eaches/box 10 boxes/case
MG100S	Medgluv Latex Exam Glove	Small	100 eaches/box 10 boxes/case
MG100M	Medgluv Latex Exam Glove	Medium	100 eaches/box 10 boxes/case
MG100L	Medgluv Latex Exam Glove	Large	100 eaches/box 10 boxes/case
MG100XL	Medgluv Latex Exam Glove	X-Large	100 eaches/box 10 boxes/case

### Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.16	6.5
Palm:	.13	5.2
Cuff:	.10	4.0
Average Length	242mm	9.5"

### Physical Properties

Before Aging	
Tensile Strength:	29 MPa
Ultimate Elongation:	700%
After Aging	
Tensile Strength:	25 MPa
Ultimate Elongation:	600%

### Quality Standards

Manufactured in accordance with Quality Sysyem ISO 9001.  
Meet and exceeds ASTM D3578 , D5151 and FDA CFR 177-2600 -  
US Food Contact Approved, QSR/FDA 510K Medical Examination Grade.

Note: Specifications are subject to change without notice.

EXCEPTIONAL FIT,  
THICKNESS FOR  
DURABLE PROTECTION



Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



# MICRO-TOUCH® NITRAFREE™

H i g h   T o u c h

## Provides comfortable and reliable hand protection

- **Reduced allergy risk:** Manufactured without sulfur-based chemical accelerators, helping to protect users from Type IV chemical allergies
- **Comfortable fit:** These nitrile medical exam gloves are made from a durable and soft nitrile formulation for enhanced comfort
- **Heightened tactility:** Thin and lightweight to ensure excellent tactility while working
- **Specialized protection:** Tested and approved for use with chemotherapy drugs\*

\*see product packaging for details

## Recommended For

- Tested for use with chemotherapy drugs
- Standard patient examinations
- Protection from latex Type I allergy and minimize chemical Type IV allergies in HCW's or patients

## FEATURED TECHNOLOGIES



## MICRO-TOUCH® Everyday, Everywhere Protection

MICRO-TOUCH® premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

# MICRO-TOUCH® NITRAFREE™

High Touch

## PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Pink	Cuff Length	Standard
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Palm Thickness (mm/mil)	0.08 / 3.1
Finger Thickness (mm/mil)	0.12 / 5.5	Allergy Prevention	Latex (Type I), Chemical (Type IV)
Available Sizes	XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Silicone Free	Not Tested	Sterile	No
Antistatic	Not Tested	Double Gloving Recommendation	Yes, as underglove
Glove Length (mm/inches)	240 / 9.5	Product Segmentation	High Touch

## ORDERING INFORMATION

Size	XS (5.5 - 6)	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034510	6034511	6034512	6034513	6034514

## PACKAGING AND STORAGE

Packaging	100 gloves per dispenser 10 dispensers per case 1000 gloves per case
Shelf Life	3 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

## STANDARDS AND CERTIFICATIONS

ASTM D6319, Category III, CE 0493, EN 16523-1, EN 421:2010, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, FDA21 CFR 177-2600-US Food Contact Approved, ISO 9001, ISO 10993-10:2002, ISO 13485

Contact your Ansell representative for ordering or more information.

North America  
US Tel: 800 952 9916  
CA Tel: 1-844-494-7854  
Email:  
insidesalesus@ansell.com

Latin America & Caribbean  
Tel: +52(442) 296 20 50  
Email: cslac@ansell.com  
Brazil Tel: +55-11-3356-3100  
Brazil Email:  
luvas.medicas@ansell.com

Europe, Middle East & Africa  
Tel: +32 (0) 2 528 74 00  
Email: info@ansell.eu

Asia Pacific  
Tel: +603 8310 6688  
Email:  
apac.medical@ansell.com

Australia & New Zealand  
Tel: +61 3 9270 7270  
Email:  
protection@ap.ansell.com

# Leather Restraint Gloves, Large



## Description

Light brown suede leather gloves have a padded lining. Black split cowhide suede glove with Kevlar lining. Gauntlet to elbow is 15"L. These gloves are NOT bite proof. Please use with caution.

Large

3M™ Avagard™ Surgical and Healthcare Hand Antiseptic with Moisturizers 16.9 fl oz, 8/cs



Product Details

Manufacturers Part Number	7000002738
Type	Gel
Capacity	16.9 oz
Package Type	Refill Bottle
Package Quantity	8
Brand	Avagard

# Antiseptic Hydrox® Topical Liquid 1 gal. Bottle



Manufacturer #	A0013
Brand	Hydrox®
Manufacturer	McKesson Medical Surgical
Country of Origin	United States
Active Ingredients	Hydrogen Peroxide
Application	Antiseptic
Container Type	Bottle
Dosage Form	Topical Liquid
Strength	3% Strength
UNSPSC Code	51473503
Volume	1 gal.



# Miltex-Integra Miltex Concentrated Liquid Ultrasonic Solution 8 oz Bottle 8oz/Bt, 12 BT/CA (3-720)



Item Code: 3-720  
Brand: Miltex-Integra Miltex

- Concentrated Liquid Ultrasonic Solution 8 oz Bottle Neutral pH Makes 32 Gallons 8oz/Bottle, 12 BT/CA

Specifications	Values
Composition	Neutral pH
Item	Ultrasonic Solution
Item Type	Concentrated Liquid
Packaging Type	Bottle
Quantity	8oz/Bottle
Secondary Quantity	Makes 32 Gallons
Volume	8 oz

## Enzymatic Instrument Cleaning Detergent, 1 Gallon



Item Code: J0649

Brand: Jorvet


### Description

This cleaner is excellent for removing stubborn organic debris that can stain your surgical instruments. It is a blend of four different enzymes to be the most effective cleaner available.

- Free of metal attacking agents. Does not contain scale forming phosphates, chlorides or hydroxides
- Concentrated for best economy. Only 1 oz. (30ml) is needed for 1 gallon of wash solution
- Cleans a wide range of metal surfaces. Safe and effective for stainless steel, aluminum, tungsten carbide, copper and brass instruments
- Available 4 gallons/case

Enzymatic Instrument Cleaning Detergent. 1 Gallon.

# IV Extension Set 30"




Item Code: J0467B

Brand: Jorvet

## Description

30" with male luer and female luer lock. 30" standard bore tubing. Sterile.

# Ketch-All Pole



Item Code: J0915B  
Brand: Jorvet

## Product Specifications

Material:	Metal, Rubber
Overall Length:	4-6'
Product Type:	Ketch All Pole
Brand:	JorVet
For Use With:	Vicious or Hard to Control Animals

# Hemastop Styptic Powder



Item Code: J1623b

Brand: Jorvet

## Description

Hemastop Styptic Powder. 1.76 oz screw cap bottle.

Effective formulation for quick stop of smaller cuts or wounds on external surfaces such as a nail trim. Use of a cotton tip applicator applied directly to wound.

Ingredients:

- p-Aminobenzenesulfonamide, Dyclonine, Plant Hemocoagulase, Desmopressin, Vitamin K, Pearl Layer Powder, Talc


## SKLAR SINGLE USE (DISPOSABLE) INSTRUMENTS

Item Description: Stitch Scissors, Econo™, Spencer-Littauer, 4 ½" Length, Floor Grade Pakistan Stainless Steel, Sterile, Finger Ring Handle, Straight Blunt Tip/Blunt Tip, Disposable, 50/cs


Supplier Item Code: 96-2596

NDC Item Code: SKR 96-2596

Unit of Measure: CS



## SYMMETRY SURGICAL SCISSORS

Item Description: Symmetry® Scissors, Operating, Straight, Sharp/Blunt, 5-1/8  
in   
Supplier Item Code: 03-0613  
NDC Item Code: AAR 03-0613  
Unit of Measure: EA



Specification sheet  
Item No:120

KIMBERLY-CLARK KIMWIPES

**KimWipes® Task Wipers, X-Large, 11.8" x 11.8", White, Disposable, Popup Box, 198/pk,**

Supplier Item Code: 34133

NDC Item Code: KIM 34133

Unit of Measure: CS





Specification sheet  
Item No:121

## Lint-Free Lens

**Item Code: J0869**

**Brand: Jorvet**

### Description

Tissues. Special grade pop-up sheets remove oil immersion that collects and dries on objectives. 3" x 4 1/2" x 4 1/2". 280/box. Similar to Kim Wipes.

## Additional information

<b>Weight</b>	<i>0.35 lbs</i>
---------------	-----------------

<b>Dimensions</b>	<i>4.75 x 4.75 x 3 in</i>
-------------------	---------------------------

<b>Microscopes &amp; Accessories</b>	<i>Cleaning</i>
--------------------------------------	-----------------



# Slide Box For 100 Slides, Cork Lined, Stainless Steel Lock, Orange

Product Description	Specifications
---------------------	----------------

Globes microscope slide storage boxes are designed for the safe long term storage of microscope slides. They are produced from durable ABS plastic and offer a hinged lid for easy opening and closing. An inventory sheet is located on the inside cover for convenient slide identification and organization. The base of the box is cork lined for added protection. Available in four sizes, these storage boxes accommodate 25 x 75mm and 1 x 3" microscope slides.

## Product Details


Manufacturers Part Number	513080N
Sterile	Non-Sterile
Color	Orange

Package Quantity	1
Brand	Globe Scientific

Specification sheet  
Item No:125

## DUKAL MVP STERILIZATION POUCHES

Item Description: Sterilization Pouches, 3.5" x 10", 200/bg, 10 bg/cs  
Supplier Item Code: 27304-BK  
NDC Item Code: DUK 27304-BK  
Unit of Measure: CS



JorVet Nylon Plastic Basket Muzzle #2



Item No: J0884A

Brand: JorVet

Description	Additional information
<h3>Description</h3> <p>This product is for muzzle use only. It is not to be used with a leash, etc to walk or restrain animal.</p> <p>Constructed of strong nylon plastic with a reinforced nose and strong leather strap. Tan.</p> <p>Size: 2 Breed: Toy, Yorkshire Terrier, Papillon, etc. End of Snout Circumference: 2" Muzzle Length: 3.5" Circumference of Bridge of Nose: 8"</p>	

Description	Additional information
<h3>Additional information</h3>	
Weight	0.1 lbs
Dimensions	5 × 2.5 × 1.5 in
Size	#2
Muzzle Type	Basket

## EXEL HYPODERMIC NEEDLES



**Item No: EXE 26417**

**Brand: Exel Corporation**

### **Description:**

Hypodermic Needle, 20G x 1", 100/bx, 20 bx/cs

#### Features/Benefits:

- Needle only
- Latex Free, This product does not contain Natural Rubber Latex
- Sterile, Non-Toxic, Non-Pyrogenic
- SUS - 304 Stainless Steel
- Needles and boxes are color-coded for easy identification
- Each needle is individually packaged, sterilized and ready for use
- Polypropylene hub will fit Luer-Lock or Luer-Slip syringes
- Hub Color: Yellow

Packages: Available 100 per box or 2000 per case

Dimensions: Length: 17 in x Width: 14 in x Height: 5 in (432 mm x 356 mm x 127 mm)

Specifications vary slightly between sizes. Packaging subject to change.

Related Categories: Needles & Syringes - Needles - Standard

#### **Needles Attributes:**

Types: Hypodermic

Configurations: Needle Only

Gauges: 20G

Lengths: 1.00 in

Specification Sheet  
Item No : 134

## **SYMMETRY SURGICAL NEEDLE HOLDER**



**Item No: AAR 10-1815**


**Brand: Symmetry Surgical, Inc.**

**Description:**

Symmetry® Needle Holder, Mayo-Hegar, 6-2/3 in, 160 mm

Specification Sheet  
Item No : 135

Cruciate Suture Needle, Large



Item No: J0710Z

Brand: Jorvet

Description	Additional information
<h3>Description</h3> <p>Cranial Cruciate Ligament Lateral Suture System. Treatment of Cranial Cruciate Ligament (CCL) rupture by placement of a lateral nylon suture has gained popularity among veterinary surgeons. The relatively quick procedure results are comparable to other techniques. Post operative strength of the repair is immediately good. Risks to intra-articular structures are minimal.</p> <p>Method: Heavy monofilament nylon is first placed on the lateral side of the stifle. A metal crimp tube replaces the large knot normally required. Appropriate tension is placed on the suture. Then a crimp tube applied over the suture ends is held in place by the crimping forceps. The lateral nylon is placed behind the lateral fabella bone and then through the tibial</p> <p>Cruciate Suture Needle. Designed for use with heavy nylon in the lateral suture repair method for cranial cruciate ruptures. Needles are specially curved along with a heavier body, making passage around the femur/fabella area much easier. Stainless steel.</p> <p>Cruciate Suture Needle. Large. 6/pk.</p>	

Description	Additional information
<h3>Additional information</h3>	
Weight	0.05 lbs
Dimensions	3.75 x 2.5 x 0.25 in
Size	Large
Orthopedic Instruments	Needles

## EXEL HYPODERMIC NEEDLES



**Item No: EXE 26411**

**Brand: EXEL**

### Description

## Hypodermic Needle 22G X 1 Inch 2000/Case 26411

EXELINT INTERNATIONAL™ is a privately held corporation that has been manufacturing and distributing disposable medical and animal health care products since the 1980's. EXELINT INTERNATIONAL's comprehensive product line features both safety and conventional disposable hypodermic products.

- 22G
- 1 Inch needle
- Hub color: Black
- This product does not contain natural rubber latex
- Sterile
- Non-toxic
- Pyrogen free
- SUS-304 stainless steel
- Each needle is individually packaged, sterilized and ready for use
- Box of 100
- Case of 20 boxes
- Made by: Exel
- Manufacturer's Product Number: 26411



Specification Sheet  
Item No : 141

Personal Eye Wash Bottle, 4 oz.



Item No: 7-006

Brand: PHYSICIANS CARE

Dimensions	
Bottle Size	4 oz
Housing Width	2 in
Housing Sz (H x W x D)	5 in x 2 in x 2 in
Housing Height	5 in
Housing Depth	2 in

Details	
Item	Single Use Eyewash Bottle
Features	Eye Flushing Solution for Irrigating the Eyes, Sterile Isotonic Buffered Solution
Shelf Life	3 yr
Standards	ANSI Z358.1-2009 (As a Personal Eye Wash Only)
Housing Included	Yes
Housing Cap.	4 oz Bottles
Activation Style	with Sterile Isotonic Buffered Solution
Brands	PHYSICIANS CARE
Housing Material	Plastic
Package Quantity	1

Is something inaccurate or missing? Help us improve our [product information](#)

Description	
Single Use Eyewash Bottle, Bottle Size 4 oz, Features Eye Flushing Solution for Irrigating the Eyes, Sterile Isotonic Buffered Solution, Housing Depth 2 in, Housing Height 5 in, Housing Width 2 in, Housing Material Plastic, Shelf Life 3 yr, Housing Included Yes, Standards ANSI Z358.1-2009 (As a Personal Eye Wash Only)	

## Jorvet Osteotome 10mm



**ITEM CODE: J0228A**

### DESCRIPTION:

Shorter length for finer control, suggested for femoral head removal. German-made. Stainless steel. 13.5cm (5").

Osteotome. Blade, 10mm W.

Description	Additional information
<b>Additional information</b>	
Weight	0.12 lbs
Dimensions	5.75 × 0.5 × 0.25 in
Length	5"
Product Material	Stainless Steel
Blade Width	10mm
Orthopedic Instruments	Osteotome

## Pet Piller, Bulk, 12/bag



**ITEM CODE: J0384**

### DESCRIPTION:


Long plastic barrel with plastic plunger. Designed like a large animal balling gun. Administer tablets and capsules to pets more quickly, safely, and easily.

- Soft latex rubber tip
- Easy to clean
- Gentle to pet's mouth
- Correct size even for small pets

Pet Piller. Bulk, 12/bag.

Description	Additional information
<b>Additional information</b>	
Weight	0.31 lbs
Dimensions	11 × 5 × 0.25 in
Pills - Pet Piller	1
Quantity	12

# Microscope Cover Slips, Plastic



ITEM CODE: J0336P

DESCRIPTION:

Almost unbreakable flexible vinyl cover slip. Interleaved with tissue. 22mm W x 22mm L x 0.157mm (0.007") H. 100/box.

Description	Additional information
-------------	------------------------

## Additional information

Weight	0.03 lbs
Dimensions	1 × 1 × 1.5 in
Laboratory Miscellaneous	Slides & Slide Holders

# Pooper Scooper



**ITEM CODE: J1103**

## DESCRIPTION:

A simple device that every clinic needs!

- Two snap together handles
- Flat side or forked, shovel side
- Bin receptacle on the second handle
- Easy to clean
- No instructions necessary

Description

Additional information

## Additional information

Weight *1.25 lbs*

Dimensions *25.75 × 8.75 × 6.5 in*

Patient Comfort *Miscellaneous*

## SIMPORT SNAPTWIST® SCINTILLATION VIALS



**ITEM CODE:** SIM S207

**Item Description:** Polyethylene Scintillation Vial, Polypropylene Cap,  
1000/cs

**Supplier Item Code:** S207

**NDC Item Code:** SIM S207

### **DESCRIPTION:**

Can be used for liquid scintillation counting, gamma counting, chromatography, sample storage and culturing. It will fit very nicely into the LKB and Packard Varisette counters. The shoulderless vial features a full width opening of 12.5mm and the overall dimensions are 16 x 57mm. The closures are made of high-flow polypropylene and are a true time saver. The vials can be securely sealed by simply snapping the caps on; removal of caps requires an easy 1/4 turn (twist). The ease with which these caps can be manipulated eliminates the danger of spillage associated with other push-on/pull-off caps. A built-in positive lock prevents the cap from popping off because of a small build-up of pressure in the vial.

## MYDENT DEFEND PROPHY PASTE



**ITEM CODE:** MYD PP-1500

**BRAND NAME:** Mydent

**Item Description:** Propy Paste, Medium-Mint, 200 cups/bx

**Supplier Item Code:** PP-1500

**NDC Item Code:** MYD PP-1500

### DESCRIPTION:

Smooth, pliable and reduced-splatter formula contains fluoride and a unique blend of cleaning and polishing agents. Safely and effectively removes stains from tooth surfaces. Gluten-free to assure against allergic reaction, DEFEND prophylaxis paste contains 1.23% fluoride ion and is available in three great tasting flavors: mint, cherry, and bubble gum. Easily rinses clean with water. 200 per package

## MYCO TECHNOCUT DISPOSABLE SCALPELS



**ITEM CODE:** MYC 6008T-15

**BRAND:** Myco Medical

**Item Description:** Scalpel, Size 15 Stainless Steel, 10/bx

**Supplier Item Code:** 6008T-15

**NDC Item Code:** MYC 6008T-15

### **DESCRIPTION:**

Individually wrapped. All plastic handles are marked with a millimeter scale. Each stainless steel blade is covered with a plastic protector. Myco TECHNOCUT® Sterile Disposable Scalpels feature blades produced from high-quality Swiss stainless steel and ergonomic, non-slip plastic handles with an integrated metric ruler. A non-slip, clear blade shield improves safety and enables easy blade size identification.



## DUKAL TECH-MED SILVER NITRATE APPLICATORS



ITEM CODE: TEC 7482

**BRAND:** DUKAL

**Item Description:** Applicator 6", 75% Silver Nitrate, 100/vl, 10 vl/bx  
(Rx) (Item is Non-Returnable) (US Only, Excluding IN  
and ND)

**Supplier Item Code:** 7482

**NDC Item Code:** TEC 7482

### DESCRIPTION:

Introducing the Silver Nitrate Applicators 6" – designed for optimal use and convenience for professionals. Our applicators are made up of 75% silver nitrate and 25% potassium nitrate, giving you top-quality materials to get the job done quickly and effectively.

The plastic caustic applicator is malleable, meaning it can be bent or shaped to access hard-to-reach areas of the ear, nose, or throat without damaging its structural integrity. And thanks to our high molded temperatures production process, every silver nitrate applicator is both hygienic and splinter-free.

Make sure your medical practice has the right tools to get any job done properly – choose Silver Nitrate Applicators 6" today!

## CONMED REFLEX® SKIN STAPLER & EXTRACTOR



ITEM CODE: CMD 8535

BRAND: Conmed Corporation

Item Description:	Skin Stapler, 35 Staples, Wide, 6/cs
Supplier Item Code:	8535
NDC Item Code:	CMD 8535

### DESCRIPTION:

The Reflex® stapler ergonomic design gives you more area between the handle and the skin (B>A). This keeps the nose of the stapler parallel to the skin allowing proper staple formation, good tissue proximation and less chance of the onset of edema.

#### **Reflex® One Skin Stapler**

- Flexibility of movement.
- Fully Retracting Anvil permits staple release while moving forward or backward
- Releases the staple when the handle is fully depressed.
- Complete staple formation is confirmed by Audible and Tactile Feedback.

## White Plastic 16 oz Slow Feeder For Cats/Dogs



ITEM CODE: 8643

BRAND: PAW PERFECT

### Description

Paw design. Helps pets eat at a slower pace. Aids in proper digestion. Helps prevent bloat. Helps prevent obesity.

### Details

Package Quantity	1
------------------	---

### Features

- Paw design & Helps pets eat at a slower pace
- Aids in proper digestion
- Helps prevent bloat
- Helps prevent obesity

PAW PERFECT

## White Plastic 16 oz Slow Feeder For Cats/Dogs



Paw design. Helps pets eat at a slower pace. Aids in proper digestion. Helps prevent bloat. Helps prevent obesity.

## HARTMANN USA ECONOLUX® GAUZE SPONGES



Woven with high quality yarn to meet low-linting, high-value needs. Ideal for use in wound cleansing, mechanical and enzyme debriding of wounds, wound packing and swabbing. 100% woven cotton and latex free. Folded edges prevent linting and unraveling. Easy opening packages.

Item                      Gauze Sponge, 2" x 2", 8-Ply, Non-Sterile, 200/slv, 25 slv/cs  
Description:  
Supplier Item    416814  
Code:

# Absorbable Hemostatic Gelatin Sponge, 10mm x 10mm



## Description

CuraSpon® Absorbable Hemostatic Gelatine Dental Cubes

- Hemostatic gelatine sponge
- Hemostasis in 2 – 4 minutes
- Liquefies in 2-5 days in contact with mucosa
- Absorbs at least 35 times its own weight in blood and fluids
- pH-neutral
- Can be used dry or impregnated
- Compatible with pharmaceuticals
- Shelf-life 4 years
- CE class III registered

CuraSpon® is manufactured from highly purified neutral gelatine foam of uniform fine porosity which guarantees a favourable hemostasis. The gelatine sponge activates the thrombocytes at the moment blood comes in contact with the matrix of the sponge. This causes the thrombocytes to release a series of substances which promote their aggregation at the same time as their surfaces change character, thus enabling them to act as a catalyst for the formation of the fibrin. CuraSpon® provides hemostasis in surgical procedures to control capillary, venous and small arterial bleeding when binding, suturing or other conventional methods of control are impractical or ineffective.

10 x 10 x 10 mm, 60 pieces



## Sterile Laparotomy Sponges 4" x 18"

10-0004

5/PK, 20 PK/CS

### Product Features

Pre-washed, 100% high quality cotton with tight weave

X-ray detectable

Banded 5s

### Attributes

Brand Name	Dukal
Size	4" x 18"
Sterile	Yes
Weight	2.64555 lbs

### Packaging

Packaging	5/PK, 20 PK/CS
Item Packaging	5/PK
Case Packaging	20 PK/CS
Total Case Quantity	100
Item Packaging Barcode	(01)10665973000081
Case Packaging Barcode	(01)20665973000088

# Feeding Dishes, Wide Rim Paw Embossed, 16 oz.



## Description

Stainless Steel:

- Heavy weight, solid, premium feeding bowls
- Heavy gauge is dishwasher safe
- Resists bacteria, stains, rust, and breakage

Wide Rim Paw Embossed Dishes. Attractive embossed paw design, light weight and less expensive.

16oz paw embossed



# Feeding Dishes, Wide Rim Paw Embossed, 32 oz.



## Description

Stainless Steel:

- Heavy weight, solid, premium feeding bowls
- Heavy gauge is dishwasher safe
- Resists bacteria, stains, rust, and breakage

Wide Rim Paw Embossed Dishes. Attractive embossed paw design, light weight and less expensive.

32oz paw embossed

AMERIHOMÉ

## Medium Stainless Steel Bucket Set, PK3



Endless uses in the house, garage, patio, and garden. Nests together for easy stacking and storage. Constructed of durable 22 gauge stainless steel. Set includes: 3 identical metal buckets: 2.37 Gallons, 9.51 Quarts, 9 Liters. 10 top, 8 bottom, 9.25 height. Long stainless handle, easy to hang and carry, 7.25 tall handle, 11.5 wide.



## Sterilization Pouches 2.25" x 4"

27301

200/BX, 20 BX/CS

### Product Features

- Self-sealing pouch has a perforated tab for easy folding and to help maintain perfect seal
- Triple sealed with extra wide "chevron" to protect instrument from piercing pouch
- Thumb notch for easy opening
- Blue see-through film allows for easy visibility inside pouch
- Rugged medical-grade paper guards against tears and perforations
- For use with steam and ETO sterilization only
- Lead free and not made with natural rubber latex

### Attributes

Brand Name	Dukal
Size	2.25" x 4"
Sterile	No
Weight	9.70034 lbs

### Packaging

Packaging	200/BX, 20 BX/CS
Inner Packaging	200/BX
Case Packaging	20 BX/CS
Total Case Quantity	4000
Item Packaging Barcode	(01)10665973027149
Inner Packaging Barcode	(01)20665973027146
Case Packaging Barcode	(01)20665973027146

## HR® LUBRICATING JELLY



HR® Lubricating Jelly is sterile, water-soluble, stable, non-greasy and spreads evenly. It is widely used in hospitals and physicians' offices...ideal for lubricating rubber products...rectal thermometers...rectal and vaginal nozzles for infants and adults. Certified Biocompatible | Micro-Flora/Sperm Friendly | Fragrance Salt and Alcohol Free | Non-Toxic/Non-Irritating & Hypoallergenic | Dye, Latex, PEG and Chlorhexidine Free | Pharmaceutical Grade Ingredients | BPOC Compliant | Medicare Reimbursement HCPCS codes: A4402 & A4332 | Duals as Ultrasound Gel | Water Soluble for quick, easy clean-up | Gluten & Glycerin Free | Kosher Certified | Made in the USA

Item	HR® Sterile Lubricating Jelly 3gm One Shot® CarePac™, 30/bx
Description:	
Supplier Item	208
Code:	

## Medline Littauer Suture Scissors, Sterile, 4-1/2"L, Stainless Steel, 50 Pack



### Features:

- Specialized scissors have a small, hook-shaped tip on one blade
- Littauer suture scissors typically have larger, less delicate tips than Spencer stitch scissors
- Stainless-steel, single-use instruments can be disposed of to help reduce risk of cross contamination

#### Product Details

Manufacturers Part Number	DYND04037	Brand	Medline
Color	Silver		

## ICU MEDICAL MEDEX STOPCOCKS



Sterile. Single. Latex free.

Item 3-Way Stopcock, Rotating Male Luer Lock,  
Description: Gamma/ETO, PVC Free, Phthalate Free, No DEHP,  
50/cs

Supplier MX5311L  
Item Code:

## TERUMO SURFLO® WINGED INFUSION SETS




The soft, flexible interlocking wings readily conform to the body's contour. The needles are bonded to the wings at just under a 3° downward angle, for less risk of vessel damage or occlusion of the lumen. The precisely honed needle assures ease of penetration and continued access with a single insertion.

Item                      Infusion Set, 21G x  $\frac{3}{4}$ ", 12" Tubing, 50/bx, 10 bx/cs (SV-  
Description:      21BLK)

Supplier Item Code:

3SV-21BLK

## TERUMO SURSHIELD® WINGED BLOOD COLLECTION & INFUSION SETS



Unique one-handed needle locking mechanism, helps reduce the risk of needle stick injuries in a cost-effective, intuitive way. The shield simply presses over the needle and locks into place with an audible click, allowing safe use without the need to change insertion technique. Sterile and disposable, adapters fit into standard infusion sets transforming them into winged blood collection sets. This allows multiple draws into evacuated tubes.

Item	Safety Winged Blood Collection Set, 21G x $\frac{3}{4}$ ", 12" Tubing, Adapter,
Description:	50/bx, 10 bx/cs (MN*SVS21B30)
Supplier Item	1MN*SVS21B30
Code:	



## TERUMO SURFLO® ETFE IV CATHETERS



An ultra-sharp double-bevel introducer needle plus a medical grade lubricant allow for easier penetration and smoother travel through tissue. Preview chamber in clear sure-grip hub gives immediate indication of proper placement. Flexible catheter with large inner diameter assures good blood flow. All this plus safe, easy-to-handle, tamper-evident packaging.

Item IV Catheter, 18G x 1¼", Green, 50/bx, 4 bx/cs (SR-  
Description: OX1832CA)

Supplier Item 3SR-OX1832CA  
Code:

# LigaFiba Lateral Suture 250lb



## Description

LigaFiba. After “fiber” an extremely strong surgical suture material. It is composed entirely of ultra high molecular weight polyethylene. It is the same hard wearing material that is used for acetubular cups in hip replacements. Weight for weight it is 15 times stronger than steel and 2.5 times stronger than nylon. It has very good abrasion characteristics which are important in hostile environments such as lateral sutures and hip toggles.

LigaFiba lateral suture. LigaFiba lateral suture with swaged on needles one end has attached cruciate needle for fabella ligament passage. The other end a straight needle for passing through tibial crest. The LigaFiba ties easily but also can be crimped (requires JorVet fiber crimper instrument J0710xa) and JorVet fiber crimp J0448 series. Sterile package.

500mm length (use J0448a crimp tube)

## MEDTRONIC POLYSORB™ COATED SYNTHETIC ABSORBABLE SUTURES



Composed of LACTOMER™ glycolide/lactide copolymer which is a synthetic polyester composed of glycolide and lactide derived from glycolic and lactic acids. POLYSORB™ sutures are prepared by coating the suture with a mixture of caprolactone/glycolide copolymer and calcium stearoyl lactylate. The advanced extrusion process of the molecule of LACTOMER™ 9-1, the exclusive braiding process and coating system give the suture: Excellent strength over the critical wound healing period; Excellent knot security; Excellent handling; Predictable absorption profile. POLYSORB™ sutures are indicated for use in soft tissue approximation or ligation and ophthalmic surgery, but not in cardiovascular or neural tissue.

Item Description: Suture, Reverse Cutting, Size 2-0, Violet, 30", Needle GS-11, 3/8 Circle, 3 dz/bx.  
Supplier Item Code: CL816

## MEDTRONIC MONOSOF™ DERMALON™ MONOFILAMENT NYLON SUTURES



MONOSOF™ monofilament nylon sutures are inert, nonabsorbable sterile surgical sutures composed of the long-chain aliphatic polymers Nylon 6 and Nylon 6.6. The advanced extrusion process of the nylon molecule gives the suture: Uniform diameter; Maximum flexibility of the strand; Excellent security with snug and flattened knots; Minimal memory; Excellent and consistent knotting strength. MONOSOF™ monofilament nylon sutures are indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological surgery.

Item	Suture, Reverse Cutting, Size 5-0, Blue, 18", Needle C-1, 3/8
Description:	Circle, 3 dz/bx
Supplier Item Code:	8886174121

## Syringe, Luer Slip, 3cc, Low Dead Space Plunger, With Cap, 100/bx, 10 bx/cs



### Features

- Sterile, Non-Toxic, Non-Pyrogenic
- Reduced dead space for reduced medication waste
- Positive plunger stops to prevent spillage when drawing
- Smooth plunger motion
- Crystal clear barrels
- Easy to read, large, bold print
- This product does not contain Natural Rubber Latex
- Available in bulk

---

### Specifications

Pieces: 1000

Weight: 14.30 lbs (6.49 kgs)

Packages: 100 per box; 1000 per case

Dimensions: Length: 28 in x Width: 12 in x Height: 11 in (711 mm x 305 mm x 279 mm)

Syringe, Luer Lock, 5-6cc, With Cap, 100/bx, 8 bx/cs (56 cs/plt)



**Item No:** EXE 26230

**Brand:** EXEL

**Product Details:**

ADDITIONAL ATTRIBUTES	Cap/Low Dead Space
DISPOSABLE/REUSABLE/REPOSABLE	Disposable
GRADUATION	Graduated
ITEM	Syringe
ITEM TYPE	General Use
LATEX/LATEX-FREE	Not Made With Natural Rubber Latex
LOW/NO DEAD SPACE	Low Dead Space
QUANTITY	100/Box
STERILE/NON-STERILE	Sterile
VOLUME	5-6cc
WITH	w/o Needle
ZONE	Nontoxic Nonpyrogenic

## General Use Syringe 3cc Low Dead Space 100/Bx, 10 BX/CA



### Specifications

Product specifications and dimensions

ADDITIONAL ATTRIBUTES	Cap
DISPOSABLE/REUSABLE/REPOSABLE	Disposable
GRADUATION	Graduated
ITEM	Syringe
ITEM TYPE	General Use
LATEX/LATEX-FREE	Not Made With Natural Rubber Latex
LOW/NO DEAD SPACE	Low Dead Space
QUANTITY	100/Box
STERILE/NON-STERILE	Sterile
VOLUME	3cc
WITH	w/o Needle
ZONE	Nontoxic Nonpyrogenic

## Syringe & Needle, Luer Lock, 3cc, Low Dead Space Plunger, 27G x 1¼", 100/bx, 10 bx/cs



Exel 3ml Syringe/Needle Combination Luer-Lock Tip is designed with a positive plunger stop to prevent spillage when drawing. The syringe is siliconized for smooth plunger action. The needle twists and locks in place on the Luer Lock Tip. Large bold print aids with greater filling accuracy. Samples may be viewed in the crystal-clear barrel.

### Features

- Sterile, Non-Toxic, Non-Pyrogenic
- Low Dead Space
- Positive plunger stops to prevent spillage when drawing
- Smooth plunger motion
- Crystal clear barrels
- Easy to read, large, bold print
- This product does not contain Natural Rubber Latex
- Available tray packaging

---

### Specifications

Pieces: 1000

Weight: 16.45 lbs (7.46 kgs)

Packages: Available 100 per box; 1000 per case

Dimensions: Length: 28 in x Width: 16 in x Height: 9 in (711 mm x 406 mm x 229 mm)



## EXEL INSULIN SYRINGE WITH NEEDLE



**Item Code: EXE 26029**

**Brand: Exel Corporation**

### Description:

Sterile. Latex free. Non-toxic. Permanently attached needle. Reduced dead space for reduced medication waste. Smooth plunger action for greater filling accuracy. Packaged in 10 packs.

Insulin Syringe & Needle, 29G x ½", 1cc, 10/bg, 10bg/bx, 5bx/cs

## Syringe, Luer Slip, 3cc, Low Dead Space Plunger, With Cap, 100/bx, 10 bx/cs



### Features

- Sterile, Non-Toxic, Non-Pyrogenic
- Reduced dead space for reduced medication waste
- Positive plunger stops to prevent spillage when drawing
- Smooth plunger motion
- Crystal clear barrels
- Easy to read, large, bold print
- This product does not contain Natural Rubber Latex
- Available in bulk

---

### Specifications

Pieces: 1000

Weight: 14.30 lbs (6.49 kgs)

Packages: 100 per box; 1000 per case

Dimensions: Length: 28 in x Width: 12 in x Height: 11 in (711 mm x 305 mm x 279 mm)

Specification Sheet  
Item No : 206

**EXEL TB TUBERCULIN SYRINGES WITH LUER SLIP**

Latex free. Sterile. 1cc syringe with low dead space plunger.

Item Description:        Tuberculin Syringe Only, 1cc, Low Dead Space  
                                 Plunger, Luer Slip with Cap, 100/bx, 10 bx/cs

Supplier Item Code:    26048



**Specifications**

Product specifications and dimensions

<b>ADDITIONAL ATTRIBUTES</b>	Cap
<b>DISPOSABLE/REUSABLE/REPOSABLE</b>	Disposable
<b>GRADUATION</b>	Graduated
<b>ITEM</b>	Syringe
<b>ITEM TYPE</b>	Tuberculin
<b>LATEX/LATEX-FREE</b>	Not Made With Natural Rubber Latex
<b>LOW/NO DEAD SPACE</b>	Low Dead Space
<b>QUANTITY</b>	100/Box
<b>STERILE/NON-STERILE</b>	Sterile
<b>VOLUME</b>	1cc
<b>WITH</b>	w/o Needle
<b>ZONE</b>	Nontoxic Nonpyrogenic

## Indicator Tape For Steam, Lead Free, .70" x 60 yds (18mm x 55m), 28/cs



Factory ISO Certification	No ISO Certification information available
Life Science Sterilization	Sterilization Indicator Tapes & Dispensers
Model	1322-18MM
Product Color	Beige
Product Type	Indicator Tape
Product Use	Exposure Monitoring
Sterilization Method	Steam
Units per Pack	28

### Dimensions and Classifications

Overall Length (Imperial)	60 yd
Overall Length (Metric)	55 m
Overall Width (Imperial)	0.709 in
Overall Width (Metric)	18 mm

**Self-Adherent Wrap, 3" x 5 yds, Tan, 24/cs**



**Item No: MMM 1583**

**Brand: 3M Health Care**

Self-adherent elastic wrap functions like a tape, but sticks only to itself. Can be used to secure, protect, compress and support. Sticks to itself; no pins or clips required. Also available with or without latex.

**Self-Adherent Wrap, 1" x 5 yds, Tan, 5/pk, 6 pk/cs**



**Item No: MMM 1581**

**Brand: 3M Health Care**

Self-adherent elastic wrap functions like a tape, but sticks only to itself. Can be used to secure, protect, compress and support. Sticks to itself; no pins or clips required. Also available with or without latex.

## Autoclave Indicator Tape, 1/2" x 60 Yd roll



### MYDENT DEFEND AUTOCLAVE INDICATOR TAPE

**Item Description:** Autoclave Indicator Tape, 1/2" x 60 Yd roll

This tape is for use in autoclave, dry heat or chemical sterilizers and will seal autoclave bags, tubing, pouches and CSR wraps. The tape will indicate when sterile conditions have been attained. Extra tacky tape to fix on all types of paper and plastic.

## Tattoo Ink, Tube Paste, Green 1oz.



### Specification:

- Non-toxic
- Tattoo Ink, Tube Paste. 1oz(140gm) Green
- Collapsible metal tube with screw cap.
- Clean area to be marked, and allow to dry
- Apply a coating of ink to the surface which is to be tattooed

### Additional information

Weight	0.05 lbs
Dimensions	5 × 1 × 1 in
Size	1 oz.
Clinic Supplies	Animal Identification



Specification Sheet  
Item No : 213

Extension Sets with MaxZero Needleless Connector



Microbore Extension. T Injection, Rem, Conn, Sic

- The MaxZero needleless connector promotes catheter patency and aids complete disinfection
  - Anti-reflux action helps prevent blood from entering the catheter upon disconnection
  - Solid, sealed surface is designed to prevent contamination
- Visible fluid path and high-flow straight fluid channels enhance flushing  
Extension sets improve convenience and flexibility

Packaging	50 Each / Case
Manufacturer	BD
Manufacturer #	MZT1003

SKU B-DMZT1003

Specifications

LATEX FREE                      Yes

## Non-Adherent Absorbent Pad, 3" x 4", Sterile, 50/bx, 18 bx/cs



### **GENTELL PRIMAPAD® NON-ADHERENT ABSORBENT PAD**

This non adherent gauze pad is often chosen as the primary dressing when extra wound bed protection is needed. Air circulates through the soft, perforated film, while fluid is absorbed by the soft cotton pad. These dressings can be easily cut to fit any shape. Does not contain natural rubber latex.



# Sterile Island Dressing 6" x 6"

4076

1/PK, 25 PK/BX, 4 BX/CS

## Product Features

- Soft, non-woven backing with adhesive border is gentle on skin conformable, breathable, and provides ease of movement
- Absorbent non-adherent pad collects excess exudates and protects wound from trauma

## Attributes

Brand Name	Dukal
Size	6" x 6"
HCPC Code	A6254
Sterile	Yes
Weight	3.17466 lbs

## Packaging

Packaging	1/PK, 25 PK/BX, 4 BX/CS
Item Packaging	1/PK
Inner Packaging	25/BX
Case Packaging	4 BX/CS
Total Case Quantity	100
Item Packaging Barcode	(01)10665973006779
Inner Packaging Barcode	(01)20665973006776
Case Packaging Barcode	(01)30665973006773

Products / Traditional Wound Care / Applicators & Swabs / Tongue Depressor / Non-Sterile Tongue Depressors 6"



## Non-Sterile Tongue Depressors 6"

9002

500/BX, 10 BX/CS

### Product Features

Clean and smooth hard wood blades

No odor or taste

Exceptional for oral use, spreading ointments, mixing medications, or splints

### Attributes

Brand Name	Dukal
Size	6"
Sterile	No
Weight	30.20333 lbs

### Packaging

Packaging	500/BX, 10 BX/CS
Inner Packaging	500/BX
Case Packaging	10 BX/CS
Total Case Quantity	5000
Inner Packaging Barcode	(01)10665973015542
Case Packaging Barcode	(01)20665973015549

## NEW WORLD IMPORTS CAREALL® TRIPLE ANTIBIOTIC



**Item ID: NWI TAO1**

**Item Description:** Triple Antibiotic Ointment, 1 oz, Compared to the Active Ingredients in Neosporin®, 24/bx

- Helps Prevent Infection in Minor Cuts, Scrapes and Burns.

**Lo-Pro Oral/ Nasal Tracheal Tube, Cuffed, Murphy Eye, 7.5mm, 10/bx**



**Item No: MAL 86052**

**Brand : Medtronic/Minimally Invasive Therapies (MIT)**

- **Pilot balloon and mechanical self-sealing valve**
- **Tip-To-Tip® radiopaque line**
- **Cuff lies close to the tube while deflated for better view of vocal cords providing ease of intubation**
- **When inflated it conforms to uneven surfaces of trachea with low intracuff pressure**
- **Low Pressure Cuff**

## PVC Stomach Tube, Small



Item No: J0106SE

Brand : Jorvet

Inexpensive standard PVC stomach tubes with funnel end.

Tube. Small. 5/16" i.d. 1/2" o.d. x 9'L

### Specification

Weight	0.65 lbs
Dimensions	9.5 × 9.5 × 0.5 in
I.D.	5/16"
O.D.	1/2"
Length	9'
Product Material	PVC
Large Animal Equipment	Stomach Tube
Size	Small

**Lo-Pro Oral/ Nasal Tracheal Tube, Cuffed, Murphy Eye, 8.0mm, 10/bx**



**Item No: MAL 86053**

**Brand: Medtronic/Minimally Invasive Therapies**

- Pilot balloon and mechanical self-sealing valve
- Tip-To-Tip® radiopaque line
- Cuff lies close to the tube while deflated for better view of vocal cords providing ease of intubation
- When inflated it conforms to uneven surfaces of trachea with low intracuff pressure



Specification Sheet  
Item No: 232

**Vacutainer® Plus Venous Blood Collection Plastic Tube, Hemogard™ Closure, 13mm x 100mm, 6.0mL, Red, Clot Activator & Silicone Coated Interior, 100/bx, 10 bx/cs**



**Item No: BEC 368044**


**Brand: BD VACUTAINER® PLUS PLASTIC BLOOD COLLECTION TUBES (SERUM)**

A safer method for blood collection. Plastic tubes reduce the risk of tube breakage and specimen spillage. Disposal of plastic tubes is safe, simple and in accordance with EPA guidelines.

## Specification

ADDITIONAL ATTRIBUTES	Silicone Coated
BRAND	Vacutainer Plus
CLOSURE	Hemogard Closure
COLOR	Red
COMPOSITION/INGREDIENTS	Plastic
DISPOSABLE/REUSABLE/REPOSABLE	Disposable
ITEM	Tube
ITEM TYPE	Venous Blood Collection
LATEX/LATEX-FREE	Not Made With Natural Rubber Latex
METHODOLOGY	Serum Chemistry

# Capillary Tube Sealant/Speciment Holding Tray



Item No: J0545

Brand: Jorvet

Tip resistant tray for upright storage of 24 tubes. Leak-tight sealant formed quickly from unique compound designed to withstand centrifugation. Non-toxic, non-drying, non-Harding.


## Specification

Weight 0.1 lbs

Dimensions 5.5 x 3.25 x 0.25 in

Tube Accessories Sealant/Holding Tray

## **Animal Health Feeding Tube, 12FR, 16", Sterile, 100/cs**



**Item No: KEN 8890701215**

**Brand: CARDINAL HEALTH**

Extra smooth uniform tubing surfaces. Smooth rounded tip.  
Thermosensitive polyvinyl chloride tubing becomes rigid when cold.

# Stomach Tube, Levin, 14FR, 48"L, Non-Sterile, 50/cs



Item Code: KEN 155711

Brand: CARDINAL HEALTH

Tube has four staggered eyes and center opening to reduce trauma. Depth marks begin 20" from distal tip with three additional marks at 4" intervals. Smooth, rounded tip reduces insertion trauma. A variety of French sizes, 48" long.

# Ultrasonic Cleaning Solution, Gallon

Item No: J0382G

Brand: Jorvet

## Description

For years, Jorgensen Labs has provided the highest quality and specialized instruments needed for the veterinary profession. Now JorVet introduces a premium line of Ultrasonic Cleaners (J0382a and J0382b) and Cleaning Solutions that keep your instruments spotless for a longer lasting performance.

What is ultrasonic cleaning and how does it work?

High frequency sound waves are used to produce more than 40,000 alternating high and low pressures per second. As these waves travel through a specially formulated cleaning solution, millions of microscopic bubbles form and violently implode. This process, known as cavitation, creates a scrubbing action capable of cleaning up to 16 times more efficiently than hand scrubbing. In minutes, all body fluids and foreign particles are removed, even from the tightest grooves, cracks and interior surfaces of the instrument being cleaned.

JorVet ultrasonic cleaning systems represent the most advanced technology in ultrasonic cleaning available today.

JorVet Concentrated Cleaning Solution is non-ammoniated. This concentrated solution is designed to dilute one part concentrate to ten parts water. It is specially formulated to provide a mild yet effective solution for cleaning instruments. Properly used, it will easily remove blood and other body fluids for longer instrument life.

JorVet Ultrasonic Cleaning Solution. Gallon.

Size

[Gallon](#)

---

## 10% Phosphate Buffered Formalin, 4 gal/cs



**Item No:** CHC C4320-101

**Brand:** Cardinal Health

**Product Details:**

- Active ingredients: formaldehyde, USP, 4%. pH is 6.9 to 7.1
- Buffered with phosphate buffer
- Not intended for use as a disinfectant

Manufacturer #	C4320-101
Brand	Cardinal Health™
Manufacturer	Cardinal
Country of Origin	TAA Compliant
Application	Histology Reagent
Chemical Name	Phosphate Buffered Formalin
Concentration Level	10%
Container Type	Bottle
Form	Solution
Type	USP Grade
UNSPSC Code	41116124
Volume	1 gal.

## SAFETY DATA SHEET

Creation Date 22-Feb-2016

Revision Date 22-Feb-2016

Revision Number 1

### 1. Identification

**Product Name** 10% Neutral Buffered Formalin

**Cat No. :** C4320-101, C4320-105, C4320-15B, C4320-180B, C4320-240B, C4320-30B, C4320-450B, C4320-60B, C4320-7.5B, C4320-90B, C4320-10H, C4320-20H, C4320-30H, C4320-45H, C4320-60H, C4322-11CCF, C4322-12CCF, C4320-10CX, C4320-20CX, C4320-30CX, C4320-30BCM, C4320-60BCM, C4320-20MD, C4320-30MD, C4320-10DL, C4320-30DL, C4320-300, C4320-10PA, C4320-20PA, C4320-30PA

**Synonyms** No information available

**Recommended Use** Laboratory chemicals.

**Uses advised against** No Information available

**Details of the supplier of the safety data sheet**

**Company**

Cardinal Health Laboratory Products  
7000 Cardinal Place Dublin, Ohio 43017  
Tel: (614) 757-7000

**Emergency Telephone Number**

Chemtrec US: (800) 424-9300  
Chemtrec EU: 001 (202) 483-7616

### 2. Hazard(s) identification

**Classification**

This chemical is considered hazardous by the 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200)

Skin Corrosion/irritation	Category 2
Serious Eye Damage/Eye Irritation	Category 1
Skin Sensitization	Category 1
Germ Cell Mutagenicity	Category 2
Carcinogenicity	Category 1A
Specific target organ toxicity (single exposure)	Category 1
Target Organs - Respiratory system, Central nervous system (CNS).	
Specific target organ toxicity - (repeated exposure)	Category 2
Target Organs - Kidney, Liver, spleen, Blood.	

**Label Elements**

**Signal Word**

Danger

**Hazard Statements**

Causes skin irritation  
May cause an allergic skin reaction

Causes serious eye damage  
May cause respiratory irritation  
May cause drowsiness or dizziness  
Suspected of causing genetic defects  
May cause cancer  
Causes damage to organs  
May cause damage to organs through prolonged or repeated exposure



### Precautionary Statements

#### Prevention

Obtain special instructions before use  
Do not handle until all safety precautions have been read and understood  
Use personal protective equipment as required  
Wash face, hands and any exposed skin thoroughly after handling  
Contaminated work clothing should not be allowed out of the workplace  
Wear protective gloves  
Do not breathe dust/fume/gas/mist/vapors/spray  
Do not eat, drink or smoke when using this product  
Use only outdoors or in a well-ventilated area

#### Response

IF exposed: Call a POISON CENTER or doctor/physician

#### Inhalation

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing  
Call a POISON CENTER or doctor/physician if you feel unwell

#### Skin

IF ON SKIN: Wash with plenty of soap and water  
Take off contaminated clothing and wash before reuse  
If skin irritation or rash occurs: Get medical advice/attention

#### Eyes

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing  
Immediately call a POISON CENTER or doctor/physician

#### Storage

Store locked up  
Store in a well-ventilated place. Keep container tightly closed

#### Disposal

Dispose of contents/container to an approved waste disposal plant

#### Hazards not otherwise classified (HNOC)

None identified

WARNING! This product contains a chemical known in the State of California to cause cancer.

### 3. Composition / information on ingredients

Component	CAS-No	Weight %
Water	7732-18-5	94 - 95
Formaldehyde	50-00-0	3.5 - 4
Methyl alcohol	67-56-1	1.2
Sodium phosphate dibasic	7558-79-4	< 1
Sodium phosphate, monobasic	7558-80-7	< 1

### 4. First-aid measures



<b>General Advice</b>	If symptoms persist, call a physician. Show this safety data sheet to the doctor in attendance.
<b>Eye Contact</b>	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Obtain medical attention. Immediately flush with plenty of water. After initial flushing, remove any contact lenses and continue flushing for at least 15 minutes. Keep eye wide open while rinsing. If symptoms persist, call a physician.
<b>Skin Contact</b>	Wash off immediately with plenty of water for at least 15 minutes. Obtain medical attention. If skin irritation persists, call a physician. Wash off immediately with soap and plenty of water while removing all contaminated clothes and shoes. If symptoms persist, call a physician.
<b>Inhalation</b>	Move to fresh air. If breathing is difficult, give oxygen. Get medical attention immediately if symptoms occur. Immediate medical attention is not required. If symptoms persist, call a physician. Move to fresh air in case of accidental inhalation of vapors or decomposition products.
<b>Ingestion</b>	Do not induce vomiting. Obtain medical attention. Clean mouth with water and drink afterwards plenty of water. Do not induce vomiting without medical advice. Never give anything by mouth to an unconscious person. Consult a physician.
<b>Most important symptoms/effects</b>	Causes eye burns. May cause allergic skin reaction. Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain or flushing
<b>Notes to Physician</b>	Treat symptomatically

## 5. Fire-fighting measures

<b>Suitable Extinguishing Media</b>	Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide. Use: Dry chemical. Carbon dioxide (CO <sub>2</sub> ). Water spray. alcohol-resistant foam.
<b>Unsuitable Extinguishing Media</b>	No information available
<b>Flash Point</b>	> 93.3 °C / > 199.9 °F
<b>Method -</b>	No information available
<b>Autoignition Temperature</b>	No information available
<b>Explosion Limits</b>	
<b>Upper</b>	No data available
<b>Lower</b>	No data available
<b>Sensitivity to Mechanical Impact</b>	No information available
<b>Sensitivity to Static Discharge</b>	No information available

### Specific Hazards Arising from the Chemical

Thermal decomposition can lead to release of irritating gases and vapors. Keep product and empty container away from heat and sources of ignition. Risk of ignition.

### Hazardous Combustion Products

Formaldehyde Methanol Carbon monoxide (CO) Carbon dioxide (CO<sub>2</sub>)

### Protective Equipment and Precautions for Firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

### NFPA

**Health**  
3

**Flammability**  
1

**Instability**  
0

**Physical hazards**  
N/A

## 6. Accidental release measures

<b>Personal Precautions</b>	Use personal protective equipment. Ensure adequate ventilation. Avoid contact with skin,
-----------------------------	--

**Environmental Precautions**

eyes and clothing. Remove all sources of ignition. Evacuate personnel to safe areas. Keep people away from and upwind of spill/leak. Pay attention to flashback. Take precautionary measures against static discharges.

Should not be released into the environment. Do not flush into surface water or sanitary sewer system. See Section 12 for additional ecological information. Prevent further leakage or spillage if safe to do so. Prevent product from entering drains.

**Methods for Containment and Clean Up** Soak up with inert absorbent material. Keep in suitable, closed containers for disposal.

## 7. Handling and storage

**Handling**

Use only under a chemical fume hood. Wear personal protective equipment. Do not breathe vapors or spray mist. Avoid contact with skin, eyes and clothing. Pay attention to flashback. Take precautionary measures against static discharges. Contents under pressure. No information available. Do not take internally. Avoid contact with clothing.

**Storage**

Keep containers tightly closed in a dry, cool and well-ventilated place. Keep containers tightly closed in a cool, well-ventilated place. Keep away from heat and sources of ignition. Keep away from heat. Keep in properly labeled containers.

## 8. Exposure controls / personal protection

**Exposure Guidelines**

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Formaldehyde	Ceiling: 0.3 ppm	(Vacated) TWA: 3 ppm (Vacated) STEL: 10 ppm (Vacated) Ceiling: 5 ppm TWA: 0.75 ppm STEL: 2 ppm	IDLH: 20 ppm TWA: 0.016 ppm Ceiling: 0.1 ppm
Methyl alcohol	TWA: 200 ppm STEL: 250 ppm Skin	(Vacated) TWA: 200 ppm (Vacated) TWA: 260 mg/m <sup>3</sup> (Vacated) STEL: 250 ppm (Vacated) STEL: 325 mg/m <sup>3</sup> Skin TWA: 200 ppm TWA: 260 mg/m <sup>3</sup>	IDLH: 6000 ppm TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> STEL: 250 ppm STEL: 325 mg/m <sup>3</sup>

Component	Quebec	Mexico OEL (TWA)	Ontario TWAEV
Formaldehyde	Ceiling: 2 ppm Ceiling: 3 mg/m <sup>3</sup>	Ceiling: 2 ppm Ceiling: 3 mg/m <sup>3</sup>	STEL: 1 ppm CEV: 1.5 ppm
Methyl alcohol	TWA: 200 ppm TWA: 262 mg/m <sup>3</sup> STEL: 250 ppm STEL: 328 mg/m <sup>3</sup> Skin	TWA: 200 ppm TWA: 260 mg/m <sup>3</sup> STEL: 250 ppm STEL: 310 mg/m <sup>3</sup>	TWA: 200 ppm STEL: 250 ppm Skin

Legend

ACGIH - American Conference of Governmental Industrial Hygienists

OSHA - Occupational Safety and Health Administration

NIOSH IDLH: The National Institute for Occupational Safety and Health Immediately Dangerous to Life or Health

**Engineering Measures**

Use only under a chemical fume hood. Ensure adequate ventilation, especially in confined areas. Ensure that eyewash stations and safety showers are close to the workstation location.

**Personal Protective Equipment****Eye/face Protection**

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166. Tightly fitting safety goggles.

<b>Skin and body protection</b>	Wear appropriate protective gloves and clothing to prevent skin exposure. Long sleeved clothing. Chemical resistant apron. Antistatic boots. Impervious gloves.
<b>Respiratory Protection</b>	Follow the OSHA respirator regulations found in 29 CFR 1910.134 or European Standard EN 149. Use a NIOSH/MSHA or European Standard EN 149 approved respirator if exposure limits are exceeded or if irritation or other symptoms are experienced.
<b>Hygiene Measures</b>	When using, do not eat, drink or smoke. Provide regular cleaning of equipment, work area and clothing.

## 9. Physical and chemical properties

<b>Physical State</b>	Liquid
<b>Appearance</b>	Clear Colorless
<b>Odor</b>	Characteristic formaldehyde
<b>Odor Threshold</b>	No information available
<b>pH</b>	7
<b>Melting Point/Range</b>	No data available
<b>Boiling Point/Range</b>	Not applicable
<b>Flash Point</b>	> 93.3 °C / > 199.9 °F
<b>Evaporation Rate</b>	No information available
<b>Flammability (solid,gas)</b>	No information available
<b>Flammability or explosive limits</b>	
Upper	No data available
Lower	No data available
<b>Vapor Pressure</b>	No information available
<b>Vapor Density</b>	No information available
<b>Specific Gravity</b>	No information available
<b>Solubility</b>	No information available
<b>Partition coefficient; n-octanol/water</b>	No data available
<b>Autoignition Temperature</b>	No information available
<b>Decomposition Temperature</b>	No information available
<b>Viscosity</b>	No information available
<b>Molecular Formula</b>	Solution

## 10. Stability and reactivity

<b>Reactive Hazard</b>	None known, based on information available
<b>Stability</b>	Stable under normal conditions.
<b>Conditions to Avoid</b>	Incompatible products. Heat, flames and sparks. Heating in air.
<b>Incompatible Materials</b>	Strong oxidizing agents, Strong acids, Strong bases
<b>Hazardous Decomposition Products</b>	Formaldehyde, Methanol, Carbon monoxide (CO), Carbon dioxide (CO <sub>2</sub> )
<b>Hazardous Polymerization</b>	Hazardous polymerization does not occur.
<b>Hazardous Reactions</b>	None under normal processing.

## 11. Toxicological information

### Acute Toxicity

<b>Product Information</b>	No acute toxicity information is available for this product
<b>Oral LD50</b>	Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.
<b>Dermal LD50</b>	Based on ATE data, the classification criteria are not met. ATE > 2000 mg/kg.
<b>Vapor LC50</b>	Based on ATE data, the classification criteria are not met. ATE > 20 mg/l.

### Component Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
-----------	-----------	-------------	-----------------

Water	-	Not listed	Not listed
Formaldehyde	500 mg/kg ( Rat )	LD50 = 270 mg/kg ( Rabbit )	0.578 mg/L ( Rat ) 4 h
Methyl alcohol	<b>Calc. ATE 60 mg/kg</b> LD50 > 1187 – 2769 mg/kg ( Rat )	<b>Calc. ATE 60 mg/kg</b> LD50 = 17100 mg/kg ( Rabbit )	<b>Calc. ATE 0.6 mg/L (vapours) or 0.5 mg/L (mists)</b> LC50 = 128.2 mg/L ( Rat ) 4 h
Sodium phosphate dibasic	LD50 = 17 g/kg ( Rat )	Not listed	Not listed
Sodium phosphate, monobasic	LD50 = 8290 mg/kg ( Rat )	LD50 > 7940 mg/kg ( Rabbit )	Not listed

**Toxicologically Synergistic Products** No information available

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

**Irritation** No information available

**Sensitization** May cause sensitization by skin contact

**Carcinogenicity** The table below indicates whether each agency has listed any ingredient as a carcinogen.

Component	CAS-No	IARC	NTP	ACGIH	OSHA	Mexico
Water	7732-18-5	Not listed	Not listed	Not listed	Not listed	Not listed
Formaldehyde	50-00-0	Group 1	Known	A2	X	A2
Methyl alcohol	67-56-1	Not listed	Not listed	Not listed	Not listed	Not listed
Sodium phosphate dibasic	7558-79-4	Not listed	Not listed	Not listed	Not listed	Not listed
Sodium phosphate, monobasic	7558-80-7	Not listed	Not listed	Not listed	Not listed	Not listed

*IARC: (International Agency for Research on Cancer)*

*IARC: (International Agency for Research on Cancer)*

*Group 1 - Carcinogenic to Humans*

*Group 2A - Probably Carcinogenic to Humans*

*Group 2B - Possibly Carcinogenic to Humans*

*NTP: (National Toxicity Program)*

*Known - Known Carcinogen*

*Reasonably Anticipated - Reasonably Anticipated to be a Human Carcinogen*

*NTP: (National Toxicity Program)*

*ACGIH: (American Conference of Governmental Industrial Hygienists)*

*A1 - Known Human Carcinogen*

*A2 - Suspected Human Carcinogen*

*A3 - Animal Carcinogen*

*ACGIH: (American Conference of Governmental Industrial Hygienists)*

*Mexico - Occupational Exposure Limits - Carcinogens*

*Mexico - Occupational Exposure Limits - Carcinogens*

*A1 - Confirmed Human Carcinogen*

*A2 - Suspected Human Carcinogen*

*A3 - Confirmed Animal Carcinogen*

*A4 - Not Classifiable as a Human Carcinogen*

*A5 - Not Suspected as a Human Carcinogen*

**Mutagenic Effects** Mutagenic effects have occurred in humans.

**Reproductive Effects** Experiments have shown reproductive toxicity effects on laboratory animals.

**Developmental Effects** Developmental effects have occurred in experimental animals.

**Teratogenicity** Teratogenic effects have occurred in experimental animals.

**STOT - single exposure** Respiratory system Central nervous system (CNS)

**STOT - repeated exposure** Kidney Liver spleen Blood

**Aspiration hazard** No information available

**Symptoms / effects, both acute and delayed** Symptoms of allergic reaction may include rash, itching, swelling, trouble breathing, tingling of the hands and feet, dizziness, lightheadedness, chest pain, muscle pain or flushing

**Endocrine Disruptor Information** No information available

**Other Adverse Effects** Tumorigenic effects have been reported in experimental animals. The toxicological

properties have not been fully investigated. See actual entry in RTECS for complete information.

## 12. Ecological information

### Ecotoxicity

Component	Freshwater Algae	Freshwater Fish	Microtox	Water Flea
Formaldehyde	Not listed	Leuciscus idus: LC50 = 15 mg/L 96h	Not listed	EC50 = 20 mg/L 96h EC50 = 2 mg/L 48h
Methyl alcohol	Not listed	Pimephales promelas: LC50 > 10000 mg/L 96h	EC50 = 39000 mg/L 25 min EC50 = 40000 mg/L 15 min EC50 = 43000 mg/L 5 min	EC50 > 10000 mg/L 24h

**Persistence and Degradability** No information available

**Bioaccumulation/ Accumulation** No information available.

### Mobility

Component	log Pow
Formaldehyde	-0.35
Methyl alcohol	-0.74

## 13. Disposal considerations

**Waste Disposal Methods** Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. Chemical waste generators must also consult local, regional, and national hazardous waste regulations to ensure complete and accurate classification.

Component	RCRA - U Series Wastes	RCRA - P Series Wastes
Formaldehyde - 50-00-0	U122	-
Methyl alcohol - 67-56-1	U154	-

## 14. Transport information

**DOT** Not regulated  
**TDG** Not regulated  
**IATA** Not regulated  
**IMDG/IMO** Not regulated

## 15. Regulatory information

**All of the components in the product are on the following Inventory lists:** Australia X = listed China Canada The product is classified and labeled according to EC directives or corresponding national laws The product is classified and labeled in accordance with Directive 1999/45/EC Europe TSCA Korea Philippines

### International Inventories

Component	TSCA	DSL	NDSL	EINECS	ELINCS	NLP	PICCS	ENCS	AICS	IECSC	KECL
Water	X	X	-	231-791-2	-		X	-	X	X	X
Formaldehyde	X	X	-	200-001-8	-		X	X	X	X	X
Methyl alcohol	X	X	-	200-659-6	-		X	X	X	X	X
Sodium phosphate dibasic	X	X	-	231-448-7	-		X	X	X	X	X
Sodium phosphate, monobasic	X	X	-	231-449-2	-		X	X	X	X	X

**Legend:**

**X** - Listed

**E** - Indicates a substance that is the subject of a Section 5(e) Consent order under TSCA.

**F** - Indicates a substance that is the subject of a Section 5(f) Rule under TSCA.

**N** - Indicates a polymeric substance containing no free-radical initiator in its inventory name but is considered to cover the designated polymer made with any free-radical initiator regardless of the amount used.

**P** - Indicates a commenced PMN substance

**R** - Indicates a substance that is the subject of a Section 6 risk management rule under TSCA.

S - Indicates a substance that is identified in a proposed or final Significant New Use Rule

T - Indicates a substance that is the subject of a Section 4 test rule under TSCA.

XU - Indicates a substance exempt from reporting under the Inventory Update Rule, i.e. Partial Updating of the TSCA Inventory Data Base Production and Site Reports (40 CFR 710(B)).

Y1 - Indicates an exempt polymer that has a number-average molecular weight of 1,000 or greater.

Y2 - Indicates an exempt polymer that is a polyester and is made only from reactants included in a specified list of low concern reactants that comprises one of the eligibility criteria for the exemption rule.

### U.S. Federal Regulations

#### TSCA 12(b)

Not applicable

#### SARA 313

Component	CAS-No	Weight %	SARA 313 - Threshold Values %
Formaldehyde	50-00-0	3.5 - 4	0.1
Methyl alcohol	67-56-1	1.2	1.0

#### SARA 311/312 Hazard Categories

Acute Health Hazard	Yes
Chronic Health Hazard	Yes
Fire Hazard	Yes
Sudden Release of Pressure Hazard	No
Reactive Hazard	No

#### CWA (Clean Water Act)

Component	CWA - Hazardous Substances	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants
Formaldehyde	X	100 lb	-	-
Sodium phosphate dibasic	X	5000 lb	-	-

#### Clean Air Act

Component	HAPS Data	Class 1 Ozone Depletors	Class 2 Ozone Depletors
Formaldehyde	X		-
Methyl alcohol	X		-

#### OSHA Occupational Safety and Health Administration

Not applicable

Component	Specifically Regulated Chemicals	Highly Hazardous Chemicals
Formaldehyde	2 ppm STEL 0.5 ppm Action Level 0.75 ppm TWA	TQ: 1000 lb

#### CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Component	Hazardous Substances RQs	CERCLA EHS RQs
Formaldehyde	100 lb	100 lb
Methyl alcohol	5000 lb	-
Sodium phosphate dibasic	5000 lb	-

#### California Proposition 65

This product contains the following proposition 65 chemicals

Component	CAS-No	California Prop. 65	Prop 65 NSRL	Category
Formaldehyde	50-00-0	Carc. (Gaseous only)	40 µg/day	Carcinogen
Methyl alcohol	67-56-1	Developmental	-	Developmental

#### U.S. State Right-to-Know Regulations

Component	Massachusetts	New Jersey	Pennsylvania	Illinois	Rhode Island
Water	-	-	X	-	-

Formaldehyde	X	X	X	X	X
Methyl alcohol	X	X	X	X	X
Sodium phosphate dibasic	X	X	X	-	-

**U.S. Department of Transportation**

Reportable Quantity (RQ): Y  
 DOT Marine Pollutant N  
 DOT Severe Marine Pollutant N

**U.S. Department of Homeland Security**

This product contains the following DHS chemicals:

Component	DHS Chemical Facility Anti-Terrorism Standard
Formaldehyde	11250 lb STQ (solution)
Sodium phosphate, monobasic	2000 lb STQ

**Other International Regulations**

Mexico - Grade Slight risk, Grade 1

**Canada**

This product has been classified in accordance with the hazard criteria of the Controlled Products Regulations (CPR) and the MSDS contains all the information required by the CPR

WHMIS Hazard Class B3 Combustible liquid  
 E Corrosive material  
 D2A Very toxic materials

**16. Other information**

Prepared By Regulatory Affairs  
 Thermo Fisher Scientific  
 Email: EMSDS.RA@thermofisher.com Richard Allan Scientific  
 A Subsidiary of Thermo Fisher Scientific  
 Tel: (800) 522-7270

Creation Date 05-Feb-2014  
 Revision Date 05-Feb-2014  
 Print Date 05-Feb-2014  
 Revision Summary This document has been updated to comply with the US OSHA HazCom 2012 Standard replacing the current legislation under 29 CFR 1910.1200 to align with the Globally Harmonized System of Classification and Labeling of Chemicals (GHS)

**Disclaimer**

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text

**End of SDS**

**Cohesive Bandage, 2" x 5 yds, Non-Sterile, Blue, Individually  
Wrapped, 36/bx**



**Item No:** AMR A2061

**Brand:** AMD-Medicom

**Product Details:**

- A highly adaptable "no-slip" bandage that sticks to itself but not to skin or other materials. Assorted colors to meet all your patient's needs. Easy one-hand tear ability. No scissors needed. Convenient dispenser box for easy storage and distribution. Medium, but adjustable, compression. Also available in blue, red, white and assorted (rainbow).



## Specification Sheet

Item No : 254

**Bandage, Cohesive, 4" x 5 yds, Latex Free (LF), Non-Sterile, Assorted Colors, 1 rl/pk, 18 pk/bx**



**Item No:** DUK 8046ASLF

**Brand:** Dukal Corporation

### **Product Details:**

Dukal 8046ASLF Cohesive Bandage

- Provides consistent support and maintains compression
- Can easily be torn by hand and is self adherent; no clips or fasteners are needed
- Lightweight, porous, and conforms well with movement
- LF series not made with natural rubber latex

Specifications of Dukal 8046ASLF Cohesive Bandage

- Brand Name Dukal
- Size 4" x 5 yd
- Color Assorted
- Contains Latex Not made with natural rubber latex
- HCPC Code A6454
- Weight 1.85 lbs