



EMPIRE

MANAGED SOLUTIONS

MEDICAL SUPPLIES Specification Catalog

Empire Medical Supplies is a subsidiary of Empire Hospitality LLC. We hold multiple GSA Contracts and are a trusted Prime Vendor for the US Military, Dept. of Veterans Affairs as well as numerous other US Government agencies. We have a proven and successful past performance with FEMA, Federal First Responder support, and state level emergency management response. We recognize the importance of ensuring the products you order meet quality and safety expectations. Our factories are established medical device manufacturers and have been producing PPE well before the COVID-19 Pandemic. We have trusted agents performing QC inspections prior to our products being shipped. Please contact us for inventory availability, delivery times, and pricing.

Sarah Jimenez

sjimenez@empiremanagementsolutions.com

615.970.2816



Contract Holder

NON-SLIP BOOT COVERS

Key Features

- Non-slip Sole
- Elastic top Opening
- Leg Ties
- Serged Seams
- Resists dry particulates, light splashes

Material Composition/Ingredients

- Polypropylene/Polyethylene
- Non-woven

Approvals and Standards

- Meets ASTM requirements

Use for:

- Splash protection.
- Disposable, Single Use

Country of Origin:

- China

Packaging:

- Carton of 50 pairs

Dimensions/Weight

- 20" height
- Thickness: 0.05mm
- Fits shoe size 8-11



Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

MEDICAL CAP



Key Features

Disposable hair cover
non-sterile
Single use
One size fits all

Material Composition/Ingredients:

Non-woven Fabric
Polypropylene SMS

Use for:

Prevent transfer of hair and skin particles.

Approvals and Standards

FDA Classification: FYF

Country of Origin:

China



Shelf life and storage

Store out of direct sunlight and away from high temperatures

Dimensions:

Package: 25cm x 25cm x 7cm
Medical Cap Circumference: 58 - 76 cm

Packaging:

Pack of 100

Special Instructions:

Dispose after use

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • (615) 970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

NITRILE GLOVES

Key Features

- Latex Free, Powder Free
- Beaded Cuff
- Ambidextrous
- Comfort fit
- Chemical Resistent



Material Composition/Ingredients

- 100% Nitrile
- Latex Free, Powder Free

Use for:

- Single Use, Disposable

Approvals and Standards

- CE Classification: Class 1 Non Sterile
- FDA Code: KHA

Physical Dimensions

Dimensions	Standards		
	Top Glove	ASTM D3578	EN 455
Length (mm)	Min 230, Min 240, 300 ± 10	Min 220 (XS, S) Min 230 (M, L, XL)	Min 240
Palm Width (mm)			
• S	84 ± 3	80 ± 10	80 ± 10
• M	94 ± 3	95 ± 10	95 ± 10
• L	105 ± 3	110 ± 10	110 ± 10
• XL	113 ± 3	120 ± 10	≥ 110
Thickness : Single Wall (mm)			
• Finger	Min 0.05	Min 0.05	N/A
• Palm	Min 0.05	Min 0.05	N/A

Shelf life and storage

Store in a cool, dry area out of direct sunlight and away from heat

Packaging:

100pc per box

Country of Origin:

China

Physical Properties

Property	ASTM D6319	EN 455
Tensile Strength (MPa)		
• Before Aging	Min 14	N/A
• After Aging	Min 14	N/A
Elongation at Break (%)		
• Before Aging	Min 500	N/A
• After Aging	Min 400	N/A
Median Force at Break (N)		
• Before Aging	N/A	Min 6
• After Aging	N/A	Min 6

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

EAR LOOP MASK



Key Features

Shingle Pleats, rectangular face masks with shapeable nose piece
Two flat knitted ear loops, one on each end

Material Composition/Ingredients

Non woven fabric
Fusible jet filter cloth
Plastic nose beam
Latex free



Approvals and Standards

FDA Class 1
Code: LYU

Shelf life and storage

Store under cool, clean, dry conditions
Avoid excessive heat (over 40°C)

Use for:

Protection from transfer of microorganisms,
& body fluids

Dimensions: Adult's Size

Length: 175 mm
Width: 90.5 mm
Ear loop: 165mm
Pleat depth: 14mm
Nose piece: 105mm

Country of Origin:

China

Dimensions: Kids' Size

Length: 145 mm
Width: 95 mm
Ear loop: 140mm
Pleat depth: 10mm
Nose piece: 95mm

Packaging:

Box of 50

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

KIDS EAR LOOP MASK

Key Features:

- Shingle pleated rectangular face mask
- Shapeable nose piece
- Two flat knitted ear loops
- Small size for young children

Material Composition/Ingredients:

- Non woven fabric
- Fusible jet filter cloth
- Plastic nose beam
- Latex Free

Approvals and Standards:

FDA Class 2

Use for:

Protect from transfer of microorganisms, & body fluids. Keep out dust, pollen and other airborne particles.

Packaging:

Individually sealed in poly bag
Box of 50

Country of Origin:

China

Dimensions/Weight:

Length: 140 mm
Width: 90 mm
Pleat depth: 10mm
Ear loops: 140 mm
Nose Piece: 105 mm

Shelf life and storage:

Store under cool, clean, dry conditions.
Avoid excessive heat (over 40 deg C)



Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

CLOTH FACE COVERING

Key Features:

Ultra soft raw edge fabric
Comfortable for extended wear
Breathable with low air-resistance
One size fits most
Easy to wash

Material Composition/Ingredients:

Thermoplastic Poly



Approvals and Standards:

Not a medically rated device

Use for:

Dust, pollen and other airborne particles

Packaging:

individually sealed in poly bag

Do not use for:

Medical Procedures/liquid transfer protection

Dimensions/Weight:

Total width: 160 mm
Ear loop hole: 63 mm
One size fits most

Country of Origin:

China

Special Instructions:

Wash after use

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

SURGICAL MASK



Key Features

Shingle Pleats, rectangular face masks with shapeable nose piece
Two flat knitted ear loops, one on each end

Material Composition/Ingredients

SMS non-woven fabric
Plastic nose beam
Latex free
Sterilized

Approvals and Standards

FFDA: Class II
Code: FXX Surgical Apparel
ASTM F2100-11 (Reapproved 2018)
ASTM F2101-14
ASTM F1862/F1862M-17
ASTM F2100-19
ASTM F2101-19

Use for:

Protection from transfer of microorganisms, body fluids & particulate material.

Country of Origin:

China



Shelf life and storage

Store under cool, clean, dry conditions
Avoid excessive heat (over 40°C)

Packaging:

Package of 50

Dimensions/Weight

Length: 12-17.5 cm
Width: 7-9.5 cm
Pleat depth: 14mm
Nose piece: 105mm

Test Method	Test Description	ASTM F2100-11 Requirement ASTM	Test Result	ASTM F2100-11 Requirement ASTM	Test Result
ASTM 1862	Fluid Resistance with synthetic blood	80 mm Hg	80 mm Hg	160 mm Hg	160 mm Hg
MIL-M-36954C	Delta P Differential Pressure, mm H ₂ O/cm ²	<4.0mm H ₂ O	3.2 mm H ₂ O	<5.0 mm H ₂ O	4.0 mm H ₂ O
ASTM F2101	Bacterial Filtration Efficiency (BFE)	≥ 95%	99.60%	≥ 98%	99.80%
ASTM F2299	Submission Particle Filtration (PFE) @ 0.1µ	≥ 95%	99.60%	≥ 98%	99.80%
16 CFR Part 1610	Flammability	Class 1	Class 1	Class 1	Class 1

Sarah Jimenez • sjimenez@empiremanagedsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

N95 MASK



Key Features

Adjustable Nose Clip

Nose Foam

Ultrasonically welded headstraps



Material Composition/Ingredients

Mask Body: SMS Fabric

Straps: Thermoplastic Elastomer

Nose Clip: Aluminum

Nose foam: Polyurethane

Contains no components made from natural rubber latex

Approvals and Standards

FDA Class 2. Code:

MSH

Use for:

Use for solid particulates and liquid mists in concentrations not exceeding 10X PEL/OEL

Do not use for:

Gases and vapors, oil aerosols, asbestos, arsenic, cadmium, lead, 4,4-methylene dianiline(MDA), or abrasive blasting

Particulate concentrations exceeding 10X PEL/OEL

Shelf life and storage

5 years from the date of manufacture

Use By date on box in MM/YYYY format

Store respirators in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals

Store in temperatures between -4°F (-20°C) and +86°F (+30°C) and not exceeding 80% RH

Dimensions/Weight

Approx weight: 0.35oz

Packaging:

Box of 10

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

N95 NIOSH MASK



Key Features

- Adjustable Nose Clip
- Nose Foam
- Ultrasonically welded head straps
- Neck worn adjustable buckle design for proper fit
- Vertical fold allows easy folding when not in use
- Fits well with goggles or helmets



Material Composition/Ingredients

- Mask Body: SMS Fabric
- Straps: Thermoplastic Elastomer
- Nose Clip: Plastic
- Nose foam: Polyurethane
- Contains no components made from natural rubber latex

Approvals and Standards

NIOSH Approved TC-84A-6660

Use for:

Use for solid particulates and liquid mists in concentrations not exceeding 10X PEL/OEL

Do not use for:

Gases and vapors, oil aerosols, asbestos, arsenic, cadmium, lead, 4,4-methylene dianiline(MDA), or abrasive blasting
Particulate concentrations exceeding 10X PEL/OEL

Shelf life and storage

- 5 years from the date of manufacture
- Use By date on box in MM/YYYY format
- Store respirators in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals
- Store in temperatures between -4°F (-20°C) and +86°F (+30°C) and not exceeding 80% RH

Dimensions/Weight

Approx weight: 0.35oz

Packaging:

Boxes of 20, 12 boxes per case

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

FACE SHIELD

Key Features

Anti Fog, Anti Static Shield
Hypoallergenic Foam head rest
Creates room for glasses or goggles
Lightweight and comfortable to wear
Splash Protection

Material Composition/Ingredients

Shield: PVC
Headrest: Foam
Headband: Elastic
Contains no natural rubber latex

Approvals and Standards

Meets ASTM requirements

Use for:

Splash protection.
Disposable, Single Use

Country of Origin:

China



Packaging:

Individually wrapped

Dimensions/Weight

Height: 9.5 inches
Width: 17.5 inches

Sarah Jimenez • sjimenez@empiremanagedsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

Key Features

Wide angle wrap around lens
Protects from sand and dust
Ventilation holes to prevent fogging
Impact resistant
Adjustable headband
Fits over most prescription glasses
Scratch resistant



Material Composition/Ingredients

Lens: High Strength Polycarbonate
Frame: PVC
Headband: Elastic

Use for:

Splash Protection

Country of Origin:

China

Approvals and Standards

ANSI Z87.1-2003
FDA Code: HOY

Packaging:

Individually Bagged
10 pcs per box

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

DISINFECTING WIPES



Key Features

- Soft damp towlette for general use
- Use on both skin and surfaces
- Resealable pack
- 75% Alcohol kills bacteria
- No Rinsing needed

Material Composition/Ingredients:

- Non-woven Fabric
- Water
- Ethanol

Use for:

Sanitizing and disinfecting surfaces at home and in public areas

Approvals and Standards

Topical Human OTC Drug
NDC: 41853-001

Country of Origin:

China

Shelf life and storage

Store out of direct sunlight and away from high temperatures

Dimensions:

Wipes: 150 mm x 200 mm
Pack: 190 mm x 110 mm x 45 mm

Packaging:

Resealable soft pack of 50 wipes.
Keep out of reach of children

Special Instructions:

Keep away from eyes and broken skin.

Sarah Jimenez • sjimenez@empiremanagedsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

HAND SANITIZER



Key Features

Clear Gel like liquid

Contains at least 70% alcohol as specified by the CDC for effective sanitizing

Material Composition/Ingredients

Principal component: 75% Dimyristyl Phosphate

Ingredient Composition Percentage

Ethanol: 70-75

Carbomer: 3-8

Glycerol: 3-8

Triethanolamine: 5-10

Water: 10-15

May contain Aloe leaf



Use for:

Skin, Hand Cleaning

Bacteriosis

Follow directions on bottle for correct usage

Shelf life and storage

Keep in properly labeled containers

Keep lid tightly closed

Store in a cool well-ventilated room

First Aid Measures:

In case of eye contact rinse thoroughly with water. Do not use on babies less than 2 months of age. Do not use on open wounds.

Discontinue use if irritation or redness occurs.

If condition persists for more than 72 hours, consult a doctor. Keep out of reach of children.

If swallowed, get medical help or contact a Poison Control Center right away.

Sizes Available:

1oz, 2oz, 4oz, 10oz, 16.9oz, 64oz

Packaging:

Clear Plastic bottles, white lids, pumps

Country of Origin:

China

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

SPRAY SANITIZER



Key Features

Antiseptic spray helps decrease bacteria
Effective at eliminating 99.9% of many common harmful germs and bacteria
No rinsing necessary

Active Ingredients:

75% v/v Ethyl Alcohol

Inactive ingredients:

Acrylates/C10-30 alkyl acrylate crosspolymer
Benzophenone-4
Fragrance
Glycerin
Tocopherol
Water

Use for:

Skin, Hand Cleaning
Sanitizing surfaces
Reducing bacteria that can cause disease

Special Instructions:

For external Use only
Keep away from eyes
May discolor fabrics

Country of Origin:

China

Approvals and Standards

Human OTC Drug
NDC 54237-001



Shelf life and storage

Keep in properly labeled containers
keep lid tightly closed
Store in a cool well-ventilated room between 68-77F
Keep away from heat, flame or sparks

Packaging:

64 oz Clear Plastic bottles with spray pump
Comes with six 10oz empty spray bottles

First Aid Measures:

In case of eye contact rinse thoroughly with water
Do not use on babies less than 2 months of age
Do not use on open wounds
Discontinue use if irritation or redness occurs
Keep out of reach of children
If swallowed seek medical advice

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

TYVEK SUIT



Key Features

- Coverall style
- Respirator-fit hood
- Elastic wrist and ankle
- Storm flap
- Zipper front closure
- Protects from dust and liquid splashes

Material Composition/Ingredients

- Tyvek, white
- Does not contain any natural rubber latex

Use for:

- Barrier protection against dry particles and aerosols
- non-hazardous light liquid splash

Do not use for:

- Contact with heavy oils, sparks or flames, or combustible liquids
- Exposure situations resulting in spray or liquid buildup on the suit
- Environments with high mechanical risks (abrasions, tears, cuts)
- Environments with exposure to hazardous substances beyond CE Type 5/6 certification
- Environments with conditions of excessive heat



Approvals and Standards

- CE approved under CE directive R 2016/425 (Personal Protective Equipment)

Sizes Available:

- S, M, L, XL, XXL

Packaging:

- Individually bagged

Country of Origin:

- China

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com

ISOLATION GOWN



Key Features

Disposable Polypropylene full length isolation Gown with ribbed SMS cuffs

Neck and waist ties

Impervious

Full Back

Autoclavable

Material Composition/Ingredients

Spun bonded Polypropylene

Comfort ribbed cuffs

Latex Free



Use for:

Single use, to protect from the transfer of microorganisms and body fluids

Shelf life and storage

Store in cool, clean dry conditions

Approvals and Standards

Class 1, Code OEA. Medical Protective Isolation Clothing. Non Sterile

Sizes Available:

S, M, L, XL, XXL

Country of Origin:

China

Packaging:

Individually bagged

Sarah Jimenez • sjimenez@empiremanagementsolutions.com • 615.970.2816

Empire Managed Solutions. 901 Officers Row, Vancouver, WA 98661. www.empireppe.com