



Upside Health Inc
Our FDA Reg: 3016685709
FDA Registered Importer & Distributor
www.upsidehealthinc.com

3-Ply Procedural Mask

ASTM Level 1 & EN14683 Type IIR Standard

Manufacturer, Importer & Product Overview

Manufacturer: Jiaxing Amazing Travel-Ware Co., Ltd.
Address: No. 2196 HONGGAO ROAD, XIUZHOU INDUSTRIAL PARK Jiaxing Zhejiang, China
FDA Registration Number: 3013888477

Spec Developer, Importer & Distributor:
Upside Health Inc
Address: 2323 Main St, Irvine CA 92614, United States of America
FDA Registration Number: 3016685709



Procedural Face Mask Testing Standards

Upside Health's 3-Ply Procedural Mask Test Results (Marked Green)						
Method	EN14683			ASTM F2100		
	Type I	Type II	Type IIR	Level 1	Level 2	Level 3
Bacterial filtration efficiency, %	≥ 95	≥ 98	≥ 98	≥ 95	≥ 98	≥ 98
Differential pressure, mm H ₂ O/cm ²	<3.0	<3.0	<5.0	<4.0	<5.0	<5.0
	<29.4	<29.4	<49.0	<39.2	<49	<49
Sub-micron particulate filtration efficiency at 0.1 micron, %	Not required	Not required	Not Required	≥ 95	≥ 98	≥ 98
Splash resistance/Synthetic Blood Resistance, mmHg Pass Result	Not required	Not required	120 (16,0 kPa)	80	120	160
Flame Spread	Not required	Not required	Not required	Class 1	Class 1	Class 1
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30	Not required	Not required	Not required



Upside Health Inc
Our FDA Reg: 3016685709
FDA Registered Importer & Distributor
www.upsidehealthinc.com

3-Ply Procedural Mask

ASTM Level 1 & EN14683 Type IIR Standard

Understanding Our Testing Standards, ASTM & EN14683

ASTM F2100 Level 1 Standards	
Characteristic	Level 1
Bacterial filtration efficiency	≥95%
Sub-micron particulates filtration efficient at 0.1 micron	≥95%
Differential pressure, mm H2O/cm2 (Breathability)	<4.0
Resistance to penetration by synthetic blood, minimum pressure in mm Hg for pass results	80 mm Hg
Flame spread	Class 1



The Importance of U.S. Based Testing: All of our ASTM F2100 Standard testing was conducted here in the United States by our trusted global partners.



International Accreditation: Our EN14683 Standard testing was conducted by internationally recognized facilities that meet Europe's gold standard.

Understanding Our Test Methods

Upside Health Protection Complies with ASTM F2100	
Specific ASTM F2100 Test Method	
ASTM F2101:	BFE ≥98%
ASTM F2299:	PFE ≥ 95% @ 0.1 micron
EN14683 C:	ΔP < 6 mmH2O/cm2
ASTM F1862	160 mmHg
16 CFR Part 1610:	Flame Spread = Class 1

Upside Health Protection Complies with EN14683 Type IIR	
Specific EN14683 Test Method	Type IIR
EN14683-5.2.2:	BFE ≥98%
EN14683 5.2.3:	ΔP < 49 Pa/cm2
EN14683-5.2.4 Splash Resistance:	16 kPa
EN14683-5.2.5 Microbial Cleanliness:	≤ 30 cfu/g



Package Configuration

SKU	Description	Qty. Bx	Qty. Cs
UPKZ-01-50	Level 1 ASTM Procedural Mask (Type IIR) – Upside Box	50	2,000
YYKZ-01-50	Level 1 ASTM Procedural Mask (Type IIR) – Jesun Box	50	2,000
YYKZ-01-20	Level 1 ASTM Procedural Mask (Type IIR) – Jesun Sleeve	20	5,000